

Product datasheet for **RC230468**

FGFR1 (NM_001174066) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR1 (NM_001174066) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	FGFR1
Synonyms:	bFGF-R-1; BFGFR; CD331; CEK; ECCL; FGFBR; FGFR-1; FLG; FLT-2; FLT2; HBGFR; HH2; HRTFDS; KAL2; N-SAM; OGD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC230468 representing NM_001174066
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTGGAGCTGGAAGTGCCTCCTCTTCTGGGCTGTGCTGGTACAGCCACTCTGCACCGCTAGGCCGT
CCCCGACCTTGCCCTGAACAAGATGCTCTCCCTCCTCGGAGGATGATGATGATGATGACTCCTCTTC
AGAGGAGAAAAGAAACAGATAACACCAAACCAAACCGTATGCCCGTAGCTCCATATTGGACATCCCCAGAA
AAGATGAAAAGAAATTGCATGCAGTGCAGGCTGCCAAGACAGTGAAGTCAAATGCCCTCCAGTGGGA
CCCCAAACCCCACTGCGCTGGTTGAAAAATGGCAAAGAATTCAAACCTGACCACAGAATTGGAGGCTA
CAAGGTCCGTTATGCCACCTGGAGCATCATAATGGACTCTGTGGTGCCTCTGACAAGGGCAACTACACC
TGCATTGTGGAGAATGAGTACGGCAGCATCAACCACACATACCAGCTGGATGTCGTGGAGCGGTCCCCTC
ACCGGCCATCCTGCAAGCAGGGTTGCCGCCAACAAAACAGTGGCCCTGGGTAGCAACGTGGAGTTCAT
GTGTAAGGTGTACAGTGACCCGACCCGACATCCAGTGGCTAAAGCACATCGAGGTGAATGGGAGCAAG
ATTGGCCAGACAACCTGCCTTATGTCCAGATCTTGAAGACTGCTGGAGTTAATACCACCGACAAAGAGA
TGGAGGTGCTTCACTTAAGAAATGTCTCCTTTGAGGACGACGGGAGTATACGTGCTTGGCGGGTAACTC
TATCGGACTCTCCATCACTCTGCATGGTTGACCGTCTTGAAGCCCTGGAAGAGAGGCCCGCAGTGATG
ACCTCGCCCTGTACCTGGAGATCATCATCTATTGCACAGGGCCCTTCTCATCTCCTGCATGGTGGGGT
CGGTATCGTCTACAAGATGAAGAGTGGTACCAAGAAGAGTGACTCCACAGCCAGATGGCTGTGCACAA
GCTGGCCAAGAGCATCCCTCTGCGCAGACAGGTAACAGTGTCTGCTGACTCCAGTGCATCCATGAACTCT
GGGTTCTTCTGGTTCGGCCATCACGGCTCTCCTCCAGTGGGACTCCCATGTAGCAGGGGTCTCTGAGT
ATGAGCTTCCCAGAACCCCTCGCTGGGAGCTGCCTCGGGACAGACTGGTCTTAGGCAAAACCCCTGGGAGA
GGGCTGCTTTGGCAGGTGGTGTGGCAGAGGCTATCGGGCTGGACAAGGACAAAACCCCAACCGTGTGACC
AAAGTGGCTGTGAAGATGTTGAAGTCCGACGCAACAGAGAAAGACTTGTGACACCTGATCTCAGAAATGG
AGATGATGAAGATGATCGGGAAGCATAAGAATATCATCAACCTGCTGGGGCCTGCACGCAGGATGGTCC
CTTGATGTATCGTGGAGTATGCCTCCAAGGGCAACCTGCGGGAGTACCTGCAGGCCCGGAGGCCCCCA
GGGCTGGAATACTGCTACAACCCAGCCACAACCCAGAGGAGCAGCTCTCCTCCAAGGACCTGGTGTCT
GCGCTACCAGGTGGCCGAGGCATGGAGTATCTGGCCTCCAAGAAGTGCATACCCGAGACCTGGCAGC
CAGGAATGTCTGGTACAGAGGACAATGTGATGAAGATAGCAGACTTTGGCCTCGCACGGGACATTAC
CACATCGACTACTATAAAAAGACAACCAACGGCCGACTGCCTGTGAAGTGGATGGCACCCGAGGCATTAT
TTGACCGGATCTACCCACCAGAGTGTGTGGTCTTTCGGGGTGTCTCTGTGGGAGATCTTCACTCT
GGGCGGCTCCCCATACCCCGGTGTGCCTGTGGAGGAACCTTTCAAGCTGCTGAAGGAGGGTACCAGCATG
GACAAGCCAGTAAGTGCACCAACGAGCTGTACATGATGATGCGGGACTGCTGGCATGCAGTGCCCTCAC
AGAGACCCACCTTCAAGCAGCTGGTGAAGACCTGGACCGCATCGTGGCCTTGACCTCCAACCAGGAGTA
CCTGGACCTGTCCATGCCCTGGACCAGTACTCCCCAGCTTCCCGACACCCGGAGCTCTACGTGCTCC
TCAGGGGAGGATCCGTCTTCTCATGAGCCGCTGCCGAGGAGCCCTGCCTGCCCCGACACCCAGCCC
AGCTTGCCAATGGCGGACTCAAACGCCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC230468 representing NM_001174066
 Red=Cloning site Green=Tags(s)

MWSWKCLLFWAVLVTATLCTARPSPTLPEQDALPSEDDDDDDSSSEKETDNTKPNRMPVAPYWTSPE
 KMEKKLHAVPAAKTVKFKCPSSGTPNPTLRWLKNGKEFKPDHRIGGYKVRATWSIIMDSVVP SDKGNYT
 CIVENEYGSINHTYQLDVVERSPHRPILQAGLPANKTVALGSNVEFMCKVYSDPQPHIQWLKHIEVNGSK
 IGPDNLPIYVQILKTAGVNTTDKEMEVLHLRNVSFEDAGEYTCLAGNSIGLSHHSAWLTVLEALEERPAM
 TSPLYLEIIYCTGAFILISCMVGSVIVYKMKSGTKKSDFFHSQMAVHKLAKSIPLRRQVTVSADSSAMNS
 GVLLVRPSRLSSSGTPMLAGVSEYELPEDPRWELPRDRLVLGKPLGEGCFQVVLAEAIGLDKDKPNRVT
 KVAVKMLKSDATEKDLSDLISEMEMMKMIGKHKNIINLLGACTQDGPLYVIVEYASKGNLREYLQARRPP
 GLEYCYNPSHNPEEQLSSKDLVSCAYQVARGMEYLASKKCIHRDLAARNVLTEDNVMKIADFLARDIH
 HIDYYKTTNGRLPVKWMapeALFDRIYTHQSDVWSFGVLLWEIFTLGGSPYPGPVVEELFKLLKEGHRM
 DKPSNCTNELYMMRDCWHAVPSQRPTFKQLVEDLDRIVALTSNQEYLDLSMPLDQYSPSPDTRSTCS
 SGEDSVFSHEPLPEEPCLPRHPAQLANGGLKRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001174066

ORF Size: 2199 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001174066.2
RefSeq Size:	4954 bp
RefSeq ORF:	2202 bp
Locus ID:	2260
UniProt ID:	P11362
Cytogenetics:	8p11.23
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Adherens junction, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton
MW:	82.2 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the fibroblast growth factor receptor (FGFR) family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds both acidic and basic fibroblast growth factors and is involved in limb induction. Mutations in this gene have been associated with Pfeiffer syndrome, Jackson-Weiss syndrome, Antley-Bixler syndrome, osteoglophonic dysplasia, and autosomal dominant Kallmann syndrome 2. Chromosomal aberrations involving this gene are associated with stem cell myeloproliferative disorder and stem cell leukemia lymphoma syndrome. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]</p>

Product images:



Circular map for RC230468