

## Product datasheet for **RC230499L3V**

### VEPH1 (NM\_001167911) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	VEPH1 (NM_001167911) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VEPH1
Synonyms:	MELT; VEPH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001167911
ORF Size:	2364 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230499).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001167911.1</a> , <a href="#">NP_001161383.1</a>
RefSeq ORF:	2367 bp
Locus ID:	79674
UniProt ID:	<a href="#">Q14D04</a>
Cytogenetics:	3q25.31-q25.32
MW:	90.1 kDa
Gene Summary:	Interacts with TGF-beta receptor type-1 (TGFBR1) and inhibits dissociation of activated SMAD2 from TGFBR1, impeding its nuclear accumulation and resulting in impaired TGF-beta signaling. May also affect FOXO, Hippo and Wnt signaling.[UniProtKB/Swiss-Prot Function]



[View online »](#)