

Product datasheet for **RC230544L3V**

CLCN2 (NM_001171088) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CLCN2 (NM_001171088) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLCN2
Synonyms:	CIC-2; cIC-2; CLC2; ECA2; ECA3; EGI3; EGI11; EGMA; EJM6; EJM8; HALD2; LKPAT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001171088
ORF Size:	2562 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230544).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001171088.2
RefSeq Size:	3167 bp
RefSeq ORF:	2565 bp
Locus ID:	1181
UniProt ID:	P51788
Cytogenetics:	3q27.1
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
MW:	93.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a voltage-gated chloride channel. The encoded protein is a transmembrane protein that maintains chloride ion homeostasis in various cells. Defects in this gene may be a cause of certain epilepsies. Four transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2012]