

Product datasheet for **RC230585L3V**

VAR2 (NM_001167733) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VAR2 (NM_001167733) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VAR2
Synonyms:	COXPD20; VALRS; VAR2L; VARSL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001167733
ORF Size:	2769 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230585).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001167733.2 , NP_001161205.1
RefSeq Size:	3394 bp
RefSeq ORF:	2772 bp
Locus ID:	57176
UniProt ID:	Q5ST30
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Aminoacyl-tRNA biosynthesis, Valine, leucine and isoleucine biosynthesis



[View online »](#)

MW: 102.4 kDa

Gene Summary: This gene encodes a mitochondrial aminoacyl-tRNA synthetase, which catalyzes the attachment of valine to tRNA(Val) for mitochondrial translation. Mutations in this gene cause combined oxidative phosphorylation deficiency-20, and are also associated with early-onset mitochondrial encephalopathies. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Aug 2014]