

Product datasheet for RC230619L3

RIM1 (RIMS1) (NM_001168407) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RIM1 (RIMS1) (NM_001168407) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	RIM1
Synonyms:	CORD7; RAB3IP2; RIM; RIM1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230619).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001168407
ORF Size:	3036 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001168407.1
RefSeq ORF:	3039 bp
Locus ID:	22999
UniProt ID:	Q86UR5
Cytogenetics:	6q13
MW:	113.4 kDa
Gene Summary:	The protein encoded by this gene is a RAS gene superfamily member that regulates synaptic vesicle exocytosis. This gene also plays a role in the regulation of voltage-gated calcium channels during neurotransmitter and insulin release. Mutations have suggested a role cognition and have been identified as the cause of cone-rod dystrophy type 7. Multiple transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Mar 2012]

Product images:



Circular map for RC230619L3