

## Product datasheet for **RC230622L2V**

### SEZ6L (NM\_001184773) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SEZ6L (NM_001184773) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SEZ6L
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001184773
ORF Size:	3072 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230622).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001184773.1</a> , <a href="#">NP_001171702.1</a>
RefSeq ORF:	3072 bp
Locus ID:	23544
UniProt ID:	<a href="#">Q9BYH1</a>
Cytogenetics:	22q12.1
Protein Families:	Druggable Genome, Transmembrane
MW:	112.1 kDa
Gene Summary:	May contribute to specialized endoplasmic reticulum functions in neurons.[UniProtKB/Swiss-Prot Function]



[View online »](#)