

Product datasheet for **RC230630L4V**

MIB2 (NM_001170687) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|----------------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | MIB2 (NM_001170687) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | MIB2 |
| Synonyms: | ZZANK1; ZZZ5 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001170687 |
| ORF Size: | 3168 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC230630). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001170687.1 , NP_001164158.1 |
| RefSeq ORF: | 2868 bp |
| Locus ID: | 142678 |
| UniProt ID: | Q96AX9 |
| Cytogenetics: | 1p36.33 |
| MW: | 114.8 kDa |
| Gene Summary: | The protein encoded by this gene is an E3 ubiquitin protein ligase that mediates ubiquitination of proteins in the Notch signaling pathway. The encoded protein may be a suppressor of melanoma invasion. [provided by RefSeq, Mar 2017] |



[View online »](#)