

Product datasheet for **RC230631L4V**

AREB6 (ZEB1) (NM_001174095) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AREB6 (ZEB1) (NM_001174095) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZEB1
Synonyms:	AREB6; BZP; DELTAEF1; FECD6; NIL2A; PPCD3; TCF8; ZFHEP; ZFHX1A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001174095
ORF Size:	3171 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230631).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001174095.1
RefSeq ORF:	3174 bp
Locus ID:	6935
UniProt ID:	P37275
Cytogenetics:	10p11.22
Protein Families:	Transcription Factors
MW:	117 kDa



[View online »](#)

Gene Summary:

This gene encodes a zinc finger transcription factor. The encoded protein likely plays a role in transcriptional repression of interleukin 2. Mutations in this gene have been associated with posterior polymorphous corneal dystrophy-3 and late-onset Fuchs endothelial corneal dystrophy. Alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2010]