

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for RC230632L3V

ABL2 (NM_001168238) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ABL2 (NM_001168238) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABL2
Synonyms:	ABLL; ARG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001168238
ORF Size:	3174 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230632).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 001168238.1</u>
RefSeq ORF:	3177 bp
Locus ID:	27
UniProt ID:	<u>P42684</u>
Cytogenetics:	1q25.2
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	ErbB signaling pathway, Viral myocarditis
MW:	116.3 kDa



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Gene Summary:This gene encodes a member of the Abelson family of nonreceptor tyrosine protein kinases.
The protein is highly similar to the c-abl oncogene 1 protein, including the tyrosine kinase,
SH2 and SH3 domains, and it plays a role in cytoskeletal rearrangements through its C-
terminal F-actin- and microtubule-binding sequences. This gene is expressed in both normal
and tumor cells, and is involved in translocation with the ets variant 6 gene in leukemia.
Multiple alternatively spliced transcript variants encoding different protein isoforms have
been found for this gene. [provided by RefSeq, Nov 2009]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US