

Product datasheet for **RC230655L4V**

SNX26 (ARHGAP33) (NM_001172630) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SNX26 (ARHGAP33) (NM_001172630) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNX26
Synonyms:	NOMA-GAP; SNX26; TCGAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001172630
ORF Size:	3369 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230655).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001172630.1
RefSeq ORF:	3372 bp
Locus ID:	115703
UniProt ID:	O14559
Cytogenetics:	19q13.12
MW:	120.2 kDa


[View online »](#)

Gene Summary:

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. Alternative splice variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Feb 2010]