

## Product datasheet for RC230665L4V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## ABL2 (NM\_001168236) Human Tagged ORF Clone Lentiviral Particle

**Product data:** 

**Product Type:** Lentiviral Particles

**Product Name:** ABL2 (NM\_001168236) Human Tagged ORF Clone Lentiviral Particle

Symbol: ABL2

Synonyms: ABLL; ARG

Mammalian Cell Puromycin

Selection:

Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM 001168236

ORF Size: 3483 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC230665).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**RefSeq:** <u>NM 001168236.1</u>

RefSeq Size: 12181 bp RefSeq ORF: 3486 bp

Locus ID: 27

 UniProt ID:
 P42684

 Cytogenetics:
 1q25.2

**Protein Families:** Druggable Genome, Protein Kinase

**Protein Pathways:** ErbB signaling pathway, Viral myocarditis





ORIGENE

**MW:** 126.7 kDa

Gene Summary: This gene

This gene encodes a member of the Abelson family of nonreceptor tyrosine protein kinases. The protein is highly similar to the c-abl oncogene 1 protein, including the tyrosine kinase, SH2 and SH3 domains, and it plays a role in cytoskeletal rearrangements through its C-terminal F-actin- and microtubule-binding sequences. This gene is expressed in both normal and tumor cells, and is involved in translocation with the ets variant 6 gene in leukemia. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, Nov 2009]