

## Product datasheet for **RC230673L4V**

### **CAMTA2 (NM\_001171166) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	CAMTA2 (NM_001171166) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CAMTA2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001171166
ORF Size:	3591 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230673).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001171166.1</a>
RefSeq ORF:	3594 bp
Locus ID:	23125
UniProt ID:	<a href="#">O94983</a>
Cytogenetics:	17p13.2
MW:	131.1 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the calmodulin-binding transcription activator protein family. Members of this family share a common domain structure that consists of a transcription activation domain, a DNA-binding domain, and a calmodulin-binding domain. The encoded protein may be a transcriptional coactivator of genes involved in cardiac growth. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jan 2010]