

Product datasheet for **RC230678L4V**

ZFAT (NM_001167583) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZFAT (NM_001167583) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZFAT
Synonyms:	AITD3; ZFAT1; ZNF406
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001167583
ORF Size:	3693 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230678).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001167583.2
RefSeq Size:	4551 bp
RefSeq ORF:	3696 bp
Locus ID:	57623
Cytogenetics:	8q24.22
MW:	137.7 kDa


[View online »](#)

Gene Summary:

This gene encodes a protein that likely binds DNA and functions as a transcriptional regulator involved in apoptosis and cell survival. This gene resides in a susceptibility locus for autoimmune thyroid disease (AITD) on chromosome 8q24. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Nov 2009]