

Product datasheet for RC230689

PARD3 (NM_001184789) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PARD3 (NM_001184789) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PARD3
Synonyms:	ASIP; Baz; PAR3; PAR3alpha; PARD-3; PARD3A; PPP1R118; SE2-5L16; SE2-5L1T1; SE2-5T2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC230689 representing NM_001184789 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAGTGACCGTGTGCTTCGGACGGACCCGGGTGGTCGTGCCGTGCGGGACGGCCACATGAAAGTTT
TCAGCCTCATCCAGCAGCGGTGACCCGCTACCGGAAGGCCATCGCCAAGGATCCAAACTACTGGATACA
GGTGCATCGCTTGAACATGGAGATGGAGGAATACTAGACCTTGATGACATTCTTTGTGATGTAGCAGAC
GATAAAGACAGACTGGTAGCAGTGTGATGAGCAGGATCCACATCACGGAGGTGATGGCACCAGTGCCA
GTTCCACGGGTACCCAGAGCCCAGAGATATTTGGTAGTGAGCTTGGCACCAACAATGTCTCAGCCTTTCA
GCCTTACCAAGCAACAAGTGAATTTGAGGTCACACCTTCAGTCTTCGAGCAAATATGCCTCTTCATGTT
CGACGCAGTAGTGACCCAGCTCTAATTGGCCTCTCCACTTCTGTGAGTGATAGTAATTTTTCTCTGAAG
AGCCTTCAAGGAAAAATCCACACGCTGGTCAACAACAGCTGGCTTCTCAAGCAGAACACTGCTGGGAG
TCCTAAAACCTGCGACAGGAAGAAAGATGAAAACACAGAAGCCTCCCGCGGGATACTAGTAAGTGGTCT
AACCAATTTAGAGAGACAATGCTCGCTCGTCTGTGAGTGCCAGTCACCCAATGGTGGGCAAGTGGCTGG
AGAAACAAGAACAGGATGAGGATGGGACAGAAGGATAACAGTCGTGTTGAACCTGTTGGACATGCTGA
CACGGGTTTGGAGCATATACCCAATTTCTCTGGATGATATGGTAAAGCTCGTAGAAGTCCCAACAGAT
GGAGGGCCTCTGGGAATCCATGTAGTGCCTTTTCTGCTCGAGGCGGCAGAACCTGGGGTTATTAGTAA
AACGATTGGAGAAAGGTGGTAAAGCTGAACATGAAAATCTTTTTCTGAGAAATGATTGCATTGTCAGGAT
TAATGATGGCGACCTTCGAAATAGAAGATTTGAACAAGCACAACATATGTTTCGCAAGCCATGCGTACA
CCCATCATTTGGTTCCATGTGGTTCTGCAGCAAATAAAGAGCAGTATGAACAACATCCCAAAGTGA
AGAACAACTACTATTCAAGCCGTTTTAGCCCTGACAGCCAGTATATTGACAACAGGAGTGTGAACAGTGC
AGGGCTTACACCGGTGCAGAGAGCACCCGACTGAACCACCCGCTGAGCAGATAGACTCTCACTCAAGA
CTACCTCATAGCGCACACCCCTCGGAAAACCACCATCCGCTCCAGCCTCGGCACCTCAGAAATGATTTA
GTACGACTGTAAGCAGTGGTTATAACACCAAAAAAATAGGCAAGAGGCTTAATATCCAGCTTAAGAAAGG
TACAGAAGGTTTGGGATTCAGCATCACTCCAGAGATGTAACAATAGGTGGCTCAGCTCCAATCTATGTG
AAAAACATTCTCCCCGGGGGGCGCCATTAGGATGGCCGACTTAAGGCAGGAGACAGACTTATAGAGG



[View online »](#)

TAAATGGAGTAGATTTAGTGGGCAAATCCCAAGAGGAAGTTGTTTCGCTGTTGAGAAGCACCAAGATGGA
AGGAACTGTGAGCCTTCTGGTCTTTCCGCCAGGAAGACGCCTTCCACCCAAGGGAAGTAAAGCAGAAGAT
GAGGATATTGTTCTTACACCTGATGGCACCAGGGAATTTCTGACATTTGAAGTCCCCTTAATGATTCAG
GATCTGCAGGCCCTTGGTGTCAAGGTAAACCGGTCAAAGAGAACCACGCAGATTTGGGAATCTT
TGCAAGTCCATTATAATGGAGGAGCAGCATCTAAAGATGGAAGGCTTCGGGTGAATGATCAACTGATA
GCAGTAAATGGAGAATCCCTGTTGGGCAAGACAAAACCAAGATGCCATGAAAACCCTAAGAAGGTCTATGT
CTACTGAAGGCAATAAACGAGGAATGATCCAGCTTATTGTTGCAAGGAGAATAAGCAAGTGAATGAGCT
GAAGTCACCTGGGAGCCCCCTGGACCTGAGCTGCCATTGAAACAGCGTTGGATGATAGAGAACGAAGA
ATTTCCATTCCCTCTACAGTGGGATTGAGGGCTTGATGAATCGCCAGCAGAAATGCTGCCCTCAGTA
GGATAATGGGTAAATACCAGCTGTCCCCTACAGTGAATATGCCCAAGATGACACTGTCATTATAGAAGA
TGACAGGTTGCCAGTGCTTCTCCACATCTCTGACCAGTCTCTTCCAGCTCCCATGATGATGTGGG
TTTGTGACGCGCAGATGCTGGTACTTGGGCAAGGCTGCAATCAGTGATTCAGCCGACTGCTCTTTGAGTC
CAGATGTTGATCCAGTCTTGTCTTTCAACGAGAAGGATTTGGACGTGAGATGCTGACGAGACTAACT
CAATACAGTGGATGACCAGAAAGCAGGTTCTCCAGCAGAGATGTGGTCTTCCCTGGGCTGAAGAAG
TCAAGCTCGTTGGAGAGTCTGCAGACCGCAGTTGCCGAGGTGACTTTGAATGGGGATATTCCTTTCCATC
GTCCACGCGCCGGATAATCAGAGGCAGGGGATGCAATGAGAGCTTCAGAGCTGCCATCGACAAATCTTA
TGATAAACCCGCGGTAGATGATGATGATGAAGGCATGGAGACCTTGGAAGAAGACACAGAAGAAAGTTCA
AGATCAGGGAGAGAGTCTGTATCCACAGCCAGTGATCAGCCTTCCCCTCTCTGGAGAGACAAATGAATG
GAAACCAAGAGAAAGGTGATAAGACTGATAGAAAAAGGATAAACTGGAAAAGAAAAGAAGAAAGATAG
AGATAAGGAGAAGGATAAAATGAAAGCCAAGAAGGGAATGCTGAAGGGCTTGGGAGACATGTTCCAGGATT
CAAGCCAAAACCGAGAATTTAGGGAACGACAAGCTCGAGAGCGTGACTATGCTGAAATTCAGATTTTC
ATCGGACATTTGGCTGTGATGATGAGTTAATGTATGGGGAGTTTCTTCTTATGAAGTTCCATGGCTCT
CAACGCTAGACCTCAGAGCCCACGAGAAGGGCATATGATGGATGCTTTGTATGCCAAGTCAAGAAGCCG
CGGAATTCAAAACCCCTACCTGTAGACAGTAACAGATCAACTCCTAGCAATCATGATCGGATACAGCGTC
TGAGGCAAGAATTTAGCAAGCAAGCAAGATGAAGATGTAGAAGATCGTCGGCGGACCTATAGTTTTGA
GCAACCTTGGCCGAACGCACGGCCGGCAGCAGAGCGGGCGACACTCGGTGTCCGTGGAGGTGCAGATG
CAGCGGCAGCGCAGGAGGAGCGCAGAGCTCCCAGCAGGCCAGCGCCAGTACAGCTCTCTGCCTCGGC
AAAGCAGGAAAAATGCCAGCTCGGTCTCCAGGACTCTTGGGAGCAGAACTACTCCCTGGGGAAGGCTT
CCAGAGTGCCAAAGAGAACCCAGTACTCCAGTACCAAGGCTCCAGGAACGGTACTCTGGGAGGACAT
GGCTTCAACGCCAGGTCATGCTGGAACTCAGGAGCTCCTTCGCCAGGAACAGAGGCGGAAGGAGCAGC
AGATGAAGAAGCAGCCTCCTTCCGAGGGGCCAGCAACTATGACTCGTATAAGAAAGTCCAGGACCCAG
TTACGCCCTCCCAAGGGGCCCTTCCGGCAAGATGTGCCCCCTCCCCTTCTCAGGTTGCGAGGCTGAAC
AGACTTCAGACTCCTGAGAAAGGAGGCCCTTCTATTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC230689 representing NM_001184789
 Red=Cloning site Green=Tags(s)

MKVTVCFGRTRVVVPCGDGHMKVFSLIQQAVTRYRKAIAKDPNYWIQVHRLEHGDGGILDLDLDDILCDVAD
 DKDRLVAVFDEQDPHHGGDGTASSTGTQSP E I F G S E L G T N N V S A F Q P Y Q A T S E I E V T P S V L R A N M P L H V
 R R S S D P A L I G L S T S V S D S N F S S E E P S R K N P T R W S T T A G F L K Q N T A G S P K T C D R K K D E N Y R S L P R D T S N W S
 N Q F Q R D N A R S S L S A S H P M V G K W L E K Q E Q D E D G T E E D N S R V E P V G H A D T G L E H I P N F S L D D M V K L V E V P N D
 G G P L G I H V V P F S A R G G R T L G L L V K R L E K G K A E H E N L F R E N D C I V R I N D G D L R N R R F E Q A Q H M F R Q A M R T
 P I I W F H V V P A A N K E Q Y E Q L S Q S E K N N Y S S R F S P D S Q Y I D N R S V N S A G L H T V Q R A P R L N H P P E Q I D S H S R
 L P H S A H P S G K P P S A P A S A P Q N V F S T T V S S G Y N T K K I G K R L N I Q L K G T E G L G F S I T S R D V T I G G S A P I Y V
 K N I L P R G A A I Q D G R L K A G D R L I E V N G V D L V G K S Q E E V S L L R S T K M E G T V S L L V F R Q E D A F H P R E L K A E D
 E D I V L T P D G T R E F L T F E V P L N D S G S A G L G V S V K G N R S K E N H A D L G I F V K S I I N G G A A S K D G R L R V N D Q L I
 A V N G E S L L G K T N Q D A M E T L R R S M S T E G N K R G M I Q L I V A R R I S K C N E L K S P G S P P G P E L P I E T A L D D R R E R
 I S H S L Y S G I E G L D E S P S R N A A L S R I M G Y Q L S P T V N M P Q D D T V I I E D D R L P V L P P H L S D Q S S S S S H D D V G
 F V T A D A G T W A K A A I S D S A D C S L S P D V D P V L A F Q R E G F G R Q I A D E T K L N T V D D Q K A G S P S R D V G P S L G L K K
 S S S L E S L Q T A V A E V T L N G D I P F H R P R P R I I R G R C N E S F R A A I D K S Y D K P A V D D D D E G M T L E E D T E E S S
 R S G R E S V S T A S D Q P S H S L E R Q M N G N Q E K G D K T D R K K D K T G K E K K K D R D K E K D M K A K K G M L K G L G D M F R I
 Q A K T R E F R E R Q A R E R D Y A E I Q D F H R T F G C D D E L M Y G G V S S Y E G S M A L N A R P Q S P R E G H M M D A L Y A Q V K K P
 R N S K P S P V D S N R S T P S N H D R I Q R L R Q E F Q Q A K Q D E D V E D R R R T Y S F E Q P W P N A R P A T Q S G R H S V S V E V Q M
 Q R Q R Q E E R E S S Q Q A Q R Q Y S S L P R Q S R K N A S S V S Q D S W E Q N Y S P G E G F Q S A K E N P R Y S S Y Q G S R N G Y L G G H
 G F N A R V M L E T Q E L L R Q E Q R R K E Q M K K Q P P S E G P S N Y D S Y K K V Q D P S Y A P P K G P F R Q D V P P S P S Q V A R L N
 R L Q T P E K G R P F Y S


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Cloning Scheme:



Plasmid Map:



ACCN: NM_001184789

ORF Size: 3819 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184789.1](#), [NP_001171718.1](#)

RefSeq ORF: 3822 bp

Locus ID: 56288

UniProt ID: [Q8TEW0](#)

Cytogenetics: 10p11.22-p11.21

Protein Pathways:	Adherens junction, Chemokine signaling pathway, Endocytosis, Neuroactive ligand-receptor interaction, Tight junction
MW:	142.2 kDa
Gene Summary:	This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]