

Product datasheet for RC230704

PARD3 (NM_001184786) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PARD3 (NM_001184786) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PARD3
Synonyms:	ASIP; Baz; PAR3; PAR3alpha; PARD-3; PARD3A; PPP1R118; SE2-5L16; SE2-5L1T1; SE2-5T2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC230704 representing NM_001184786 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGAAAGTGACCGTGTGCTTCGGACGGACCCGGGTGGTCGTGCCGTGCGGGGACGGCCACATGAAAGTTT
TCAGCCTCATCCAGCAGGCGGTGACCCGCTACCGGAAGGCCATCGCCAAGGATCCAACTACTGGATACA
GGTGCATCGTTGGAACATGGAGATGGAGGAATACTAGACCTTGATGACATTCTTTGTGATGTAGCAGAC
GATAAAGACAGACTGGTAGCAGTGTGGATGAGCAGGATCCACATCACGGAGGTGATGGCACCAGTGCCA
GTTCCACGGGTACCCAGAGCCCAGAGATATTTGGTAGTGAGCTTGGCACCAACAATGTCTCAGCCTTTCA
GCCTTACCAAGCAACAAGTAAAATTGAGGTCACACCTTCAGTCCTTCGAGCAAATATGCCTTTCATGTT
CGACGCAGTAGTGACCCAGCTCTAATTGGCCTCTCCACTTCTGTGAGTATAGTAATTTTTCTCTGAAG
AGCCTTCAAGGAAAAATCCCACACGCTGGTCAACAACAGCTGGCTTCTCAAGCAGAACACTGCTGGGAG
TCCTAAAACCTGCGACAGGAAGAAAGATGAAAACACAGAACCTCCCGCGGATACTAGTAACCTGGTCT
AACCAATTTAGAGAGACAATGCTCGCTCGTCTGAGTGCCAGTCACCAATGGTGGGCAAGTGGCTGG
AGAAACAAGAACAGGATGAGGATGGGACAGAAGAGGATAACAGTCGTGTTGAACCTGTTGGACATGCTGA
CACGGGTTTGGAGCATATACCAACTTTTCTCTGGATGATATGGTAAAGCTCGTAGAAGTCCCAACGAT
GGAGGGCTCTGGGAATCCATGTAGTGCCTTTTCAGTGCTCGAGGGCGCAGAACCTGGGGTTATTAGTAA
AACGATTGGAGAAAGGTGTAAGCTGAACATGAAAATCTTTTTCTGTGAGAATGATTGCATTGTCAGGAT
TAATGATGGCGACCTTCGAAATAGAAGATTTGAACAAGCACAAACATATGTTTCGCAAGCCATGCGTACA
CCCATCATTTGGTTCCATGTGGTTCCTGCAGCAAATAAGAGCAGTATGAACAACATCCCAAAGTGA
AGAACAATTAATAAAGCCGTTTTAGCCCTGACAGCCAGTATATTGACAACAGGAGTGTGAACAGTGC
AGGGCTTACACGGTGCAGAGAGCACCCCGACTGAACCACCCGCTGAGCAGATAGACTCTCACTCAAGA
CTACCTCATAGCGCACACCCCTCGGAAAACCACCATCCGCTCCAGCCTCGGCACCTCAGAATGTATTTA
GTACGACTGTAAGCAGTGGTTATAACACCAAAAAAATAGGCAAGAGGCTTAATATCCAGCTTAAGAAAGG
TACAGAAGGTTTGGGATTCAGCATCACTCCAGAGATGTAACAATAGGTGCTCAGTCCAATCTATGTG



[View online >](#)

AAAAACATTCTCCCCGGGGGGCGGCCATTTCAGGATGGCCGACTTAAGGCAGGAGACAGACTTATAGAGG
TAAATGGAGTAGATTTAGTGGGCAAATCCCAAGAGGAAGTTGTTTCGCTGTTGAGAAGCACCAAGATGGA
AGGAACTGTGAGCCTTCTGGTCTTTTCGCCAGGAAGACGCCTTCCACCCAAGGAACTGAAAGCAGAAGAT
GAGGATATTGTTCTTACACCTGATGGCACCAGGGAATTTCTGACATTTGAAGTCCCCTTAATGATTCAG
GATCTGCAGGCCCTGGTGTCAAGGTAAACCGGTCAAAGAGAACCACGCAGATTTGGGAATCTT
TGTC AAGTCCATTATTAATGGAGGAGCAGCATCTAAAGATGGAAGGCTTCGGGTGAATGATCAACTGATA
GCAGTAAATGGAGAATCCCTGTTGGCAAGCAAACCAAGATGCCATGGAACCCCTAAGAAGGTCTATGT
CTACTGAAGGCAATAAACGAGGAATGATCCAGCTTATTGTTGCAAGGAGAATAAGCAAGTCAATGAGCT
GAAGTCACCTGGGAGCCCCCTGGACCTGAGCTGCCATTGAAACAGCGTTGGATGATAGAGAACGAAGA
ATTTCCCATTCCTCTACAGTGGGATTGAGGGCTTGATGAATCGCCAGCAGAAATGCTGCCCTCAGTA
GGATAATGGGTAAATACCAGCTGTCCCCTACAGTGAATATGCCCAAGATGACTGTCTTATAGAAGA
TGACAGGTTGCCAGTGCTTCTCCACATCTCTGACCAGTCTCTTCCAGCTCCCATGATGATGTGGG
TTTGTGACGCAGATGCTGGTACTTGGCCAAGGCTGCAATCAGTGATTCAGCCGACTGCTTTGAGTC
CAGATGTTGATCCAGTCTTGTCTTTCAACGAGAAGGATTTGGACGTGAGATGTCAGAAAAACGCAC
AAAGCAATTTTCAGATGCCAGTCAATTGGATTTGTTAAAACACGAAAATCAAAAAGCATGGATTTAGGT
ATAGCTGACGAGACTAACTCAATACAGTGGATGACCAGAAAGCAGGTTCTCCACAGCAGATGTGGGTC
CTTCCCTGGGTCTGAAGAAGTCAAGCTCGTTGGAGAGTCTGCAGACCGCAGTTGCCGAGGTGACTTTGAA
TGGGGATATTCCTTTCCATCGTCCACGGCCCGGATAATCAGAGGCAGGGGATGCAATGAGAGCTTCAGA
GCTGCCATCGACAAATCTTATGATAAACCCGCGGTAGATGATGATGAAGGCATGGAGACCTTGGAA
AAGACACAGAAGAAAGTTCAAGATCAGGGAGAGAGTCTGTATCCACAGCCAGTGATCAGCCTTCCCCTC
TCTGGAGAGACAAATGAATGGAAACCAAGAGAAAGGTGATAAGACTGATAGAAAAAGGATAAACTGGA
AAAGAAAAGAAGAAAGATAGAGATAAGGAGAAGGATAAAATGAAAGCCAAGAAGGGAATGCTGAAGGGCT
TGGGAGACATGTTTCAGGTTTGGCAAACATCGAAAAGATGACAAGATTGAGAAAACGGGTAATAAAAAAT
ACAGGAATCCTTTACATCAGAAGAGGAGAGGATACGAATGAAGCAGGAGCAGGAGAGGATTCAAGCCAAA
ACTCGAGAATTTAGGGAACGACAAGCTCGAGAGCGTGACTATGCTGAAATTCAGATTTTCATCGGACAT
TTGGCTGTGATGATGAGTTAATGTATGGGGGAGTTTCTTCTATGAAGGTTCCATGGCTCTCAACGCTAG
ACCTCAGAGCCCACGAGAAGGGCATATGATGGATGCTTTGTATGCCCAAGTCAAGAAGCCGCGGAATTCC
AAACCCTCACCTGTAGACAGTAACAGATCAACTCCTAGCAATCATGATCGGATACAGCGTCTGAGGCAAG
AATTTTCAGCAAGCAAAGCAAGATGAAGATGTAGAAGATCGTCGGCGGACCTATAGTTTTGAGCAACCCTG
GCCGACGCACGGCCGGCGACGCAGAGCGGGCGACTCGGTGTCCGTGGAGGTGCAGATGCAGCGGCAG
CGGCAGGAGGAGCGGAGAGCTCCACGAGGCCAGCGCCAGTACAGCTCTCTGCCTCGGCAAAGCAGGA
AAAATGCCAGCTCGGTCTCCAGGACTCTTGGGAGCAGAACTACTCCCCTGGGGAAGGCTTCCAGAGTGC
CAAAGAGAACCCAGGTAATCCAGCTACCAAGGCTCCAGGAACGGCTACCTGGGAGGACATGGCTTCAAC
GCCAGGTCATGCTGAAACTCAGGAGCTCTTTCGCCAGGAACAGAGGCGGAAGGAGCAGCAGATGAAGA
AGCAGCTCCTTCCGAGGGGCCAGCAACTATGACTCGTATAAGAAAGTCCAGGACCCAGTTACGCCCC
TCCAAGGGGCCCTTCCGGCAAGATGTGCCCCCTCCCCTTCTCAGGTTGCGAGGCTGAACAGACTTCAG
ACTCCTGAGAAAGGGAGGCCCTTCTATTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC230704 representing NM_001184786
 Red=Cloning site Green=Tags(s)

```

MKVTVCFGRTRVVVPCGDGHMKVFSLIQQAVTRYKAIKDPNYWIQVHRLHEHGGGILDLDLDDILCDVAD
DKDRLVAVFDEQDPHHGGDGTASSTGTQSP E I F G S E L G T N N V S A F Q P Y Q A T S E I E V T P S V L R A N M P L H V
R R S S D P A L I G L S T S V S D S N F S S E E P S R K N P T R W S T T A G F L K Q N T A G S P K T C D R K K D E N Y R S L P R D T S N W S
N Q F Q R D N A R S S L S A S H P M V G K W L E K Q E Q D E D G T E E D N S R V E P V G H A D T G L E H I P N F S L D D M V K L V E V P N D
G G P L G I H V V P F S A R G G R T L G L L V K R L E K G G K A E H E N L F R E N D C I V R I N D G D L R N R R F E Q A Q H M F R Q A M R T
P I I W F H V V P A A N K E Q Y E Q L S Q S E K N N Y S S R F S P D S Q Y I D N R S V N S A G L H T V Q R A P R L N H P P E Q I D S H S R
L P H S A H P S G K P P S A P A S A P Q N V F S T T V S S G Y N T K K I G K R L N I Q L K G T E G L G F S I T S R D V T I G G S A P I Y V
K N I L P R G A A I Q D G R L K A G D R L I E V N G V D L V G K S Q E E V S L L R S T K M E G T V S L L V F R Q E D A F H P R E L K A E D
E D I V L T P D G T R E F L T F E V P L N D S G S A G L G V S V K G N R S K E N H A D L G I F V K S I I N G G A A S K D G R L R V N D Q L I
A V N G E S L L G K T N Q D A M E T L R R S M S T E G N K R G M I Q L I V A R R I S K C N E L K S P G S P P G P E L P I E T A L D D R R E R
I S H S L Y S G I E G L D E S P S R N A A L S R I M G Y Q L S P T V N M P Q D D T V I I E D D R L P V L P P H L S D Q S S S S S H D D V G
F V T A D A G T W A K A A I S D S A D C S L S P D V D P V L A F Q R E G F G R Q S M S E K R T K Q F S D A S Q L D F V K T R K S K S M D L G
I A D E T K L N T V D D Q K A G S P S R D V G P S L G L K K S S S L E S L Q T A V A E V T L N G D I P F H R P R P R I I R G R C N E S F R
A A I D K S Y D K P A V D D D D E G M E T L E E D T E E S S R S G R E S V S T A S D Q P S H S L E R Q M N G N Q E K G D K T D R K K D K T G
K E K K K D R D K E K D K M K A K K G M L K G L G D M F R F G K H R K D D K I E K T G K I K I Q E S F T S E E E R I R M K Q E Q E R I Q A K
T R E F R E R Q A R E R D Y A E I Q D F H R T F G C D D E L M Y G V S S Y E G S M A L N A R P Q S P R E G H M M D A L Y A Q V K K P R N S
K P S P V D S N R S T P S N H D R I Q R L R Q E F Q Q A K Q D E D V E D R R R T Y S F E Q P W P N A R P A T Q S G R H S V S V E V Q M Q R Q
R Q E E R E S S Q Q A Q R Y S S L P R Q S R K N A S S V S Q D S W E Q N Y S P G E G F Q S A K E N P R Y S S Y Q G S R N G Y L G G H G F N
A R V M L E T Q E L L R Q E Q R R K E Q Q M K K Q P P S E G P S N Y D S Y K K V Q D P S Y A P P K G P F R Q D V P P S P S Q V A R L N R L Q
T P E K G R P F Y S
  
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja3369_c02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




* The last codon before the Stop codon of the ORF

ACCN: NM_001184786

ORF Size: 4020 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001184786.1 , NP_001171715.1
RefSeq Size:	5965 bp
RefSeq ORF:	4023 bp
Locus ID:	56288
UniProt ID:	Q8TEW0
Cytogenetics:	10p11.22-p11.21
Protein Pathways:	Adherens junction, Chemokine signaling pathway, Endocytosis, Neuroactive ligand-receptor interaction, Tight junction
MW:	149.7 kDa
Gene Summary:	This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]

Product images:



Circular map for RC230704