

## Product datasheet for **RC230725**

### **RNF17 (NM\_001184993) Human Tagged ORF Clone**

#### Product data:

|                           |                                                                                              |
|---------------------------|----------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                          |
| Product Name:             | RNF17 (NM_001184993) Human Tagged ORF Clone                                                  |
| Tag:                      | Myc-DDK                                                                                      |
| Symbol:                   | RNF17                                                                                        |
| Synonyms:                 | Mmip-2; SPATA23; TDRD4                                                                       |
| Mammalian Cell Selection: | Neomycin                                                                                     |
| Vector:                   | pCMV6-Entry (PS100001)                                                                       |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                         |
| ORF Nucleotide Sequence:  | <p>&gt;RC230725 representing NM_001184993</p> <p>Red=Cloning site Blue=ORF Green=Tags(s)</p> |

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGC**C

ATGGCGGCAGAGGCTTCGAAGACTGGGCCTTCTAGGTCTTCTACCAGCGAATGGGGAGGAAGAGTCAGC  
 CCTGGGGTGCCGCTGAAATCCAGTGCACCAGGTGTGGAAGGAGGGTATCCAGATCATCCGGTCACCATTTG  
 TGAACCTCAATGTGGACATGCTTTTGTGAACATGCTTGTTAATGACTGAAGAATGCACCACAATTATA  
 TGCCCTGATTGTGAGGTGCTACAGCTGTAATACTAGACAACGCTACTACCAATGGCTGGATATATTA  
 AGGAAGACTCCATAATGGAAAACTGCAGCCTAAGACGATAAAGAATTGTTCTCAGGACTTTAAGAAGAC  
 TGCTGATCAGCTAACTACTGGTTTGAACGTTACGCCTCCACAGACAAGACTCTTTGAACATCATCAGCT  
 GTAATGTTGGACACTAATACTGCAGAAGAAATTGATGAAGCATTGAATACAGCACACCATAGTTTGAAC  
 AGTTAAGCATTGCTGGAAAAGCACTTGAACACATGCAGAAGCAAACGATAGAGGAAAGAGAAAGAGTTAT  
 AGAAGTTGTGGAGAAACAGTTTGACCAACTTTTGGCTTTTTTGTATCCAGGAAAAAGAACCTGTGTGAA  
 GAATTTGCAAGAACTACTGATGATTATCTATCAAATTTAATAAAGGCTAAAGCTACATTGAAGAGAAAA  
 AAAATAATTTGAATGCAGCTATGAACATAGCAAGAGCATTACAATTATCGCCTTCTCTAAGAACATACTG  
 TGACCTGAATCAGATTATCCGGACTTTGCAGTTAACTTCAGATAGTGAATTAGCACAAGTTAGTTCTCCA  
 CAACTAAGGAACCCTCCAGGTTGAGTGTGAATTGCAGTGAGATCATCTGTATGTTCAACAATATGGGAA  
 AGATTGAATTTAGGGACTCAACAAAATGTTATCCCAAGAAAATGAAATTAGACAGAATGTTCAAAAAGAA  
 ATATAATAACAAAAGGAACTTTCTTGTACGATACATACCCACCGCTAGAAAAGAAAAAGGTTGACATG  
 TCTGTCTAACCAGTGAAGCACCACCACCTCCTTTGCAACCTGAGACAAATGATGTACATTTAGAAGCAA  
 AAAACTTCCAGCCACAGAAAGACGTTGCAACAGCATCCCTAAAACCATGCTGTGTACCTCAGATGGG  
 ATCTAGCCCTGATGTGATAATTGAAGAAATTATTGAAGACAACGTGAAAGTTCTGCAGAGCTAGTTTTT  
 GTAAGCCATGTAATAGATCCTTGCCATTTCTACATTGGAAGTATTCACAAATAAAGACGCCAAAGTAC  
 TGGAGAAGAAGGTGAATGAATTTTGAATAGGAGTTCACACCTTGATCCTTCAGACATTTTGAAGTACG  
 TGCAAGAATATTTGTCAGCAGTATTAATAATGGAATGTGGTGTGAGGAAGTATCACAGAATTAATTCCA


[View online »](#)

ATAGAGGGTAGAAATACCAGAAAACCTTGTAGTCCAACCAGATTATTTGTCCATGAAGTTGCACTAATAC  
 AAATATTCATGGTAGATTTTGGAAATCTGAAGTCCTGATTGTCACTGGAGTTGTTGATACCCATGTGAG  
 ACCAGAACTCTGCTAAGCAACATATTGCACTAAATGATTTATGTCTGGTTCTAAGGAAATCTGAACCA  
 TATACTGAAGGGCTGCTAAAAGACATCCAGCCATTAGCACAAACCATGCTCATTGAAAGACATTGTTCCAC  
 AGAATTCAAATGAAGGCTGGGAAGAGGAAGCTAAAGTGGAATTTTGGAAATGGTAAATAACAAGGCTGT  
 TTCAATGAAAGTTTTAGAGAAGAAGATGGTGTGCTTATTGTAGATCTGCAAAAACCCACCCTGAATAAA  
 ATAAGCAGTGATATGCCTGTGTCTCTTAGAGATGCGCTAGTTTTATGGAAGTACGAAAGTTAAGTCAC  
 AATCACTAAGAAGTCACCTTTGAAAAAAATACTACTTTACACTATCATCCACCTATTTTGCCTAAAGAAAT  
 GACAGATGTTTCAGTAACGGTTTGTATATAAATAGTCCTGGAGATTTCTATCTTCAGTTGATAGAGGGC  
 CTGGATATTTTATTTCTATTAAGACAATCGAGGAATTCTATAAAGTGAAGATGGAGAAAATCTGGAAA  
 TCCTCTGTCCAGTTCAAGATCAAGCCTGTGTAGCTAAATTTGAAGATGGAATTTGGTACCGAGCAAAAGT  
 TATCGGATTGCCTGGACATCAGGAAGTTGAAGTTAAATATGTGGACTTTGGTAATACTGCAAAAATAACA  
 ATCAAAGACGTGCGTAAAATAAGGATGAGTTTCTGAATGCCCCAGAGAAGGCAATTAATGTAAGTTGG  
 CCTATATTGAACCATATAAAGGACAATGCAGTGGTCCAAAGAAGCTAAAGAAAAATTTGAAGAAAAGGC  
 TCAAGATAAAATTTATGACATGTTCAAGTTATCAAAATCTGGAAGATAATGTGCTCTTAGTTGAGCTTTTC  
 GATTCTCTTGGTGCTCCTGAAATGACTACTACTAGTATTAATGACCAGCTAGTTAAAGAGGGCCTAGCAT  
 CTTATGAAATAGGATACATCCTCAAAGATAATTCTCAAAGCATATTGAAGTTTGGGATCCTTCTCCAGA  
 AGAAATATTTCAAATGAAGTACACAACCTAAATCCTGTGTCTGCAAAATCTCTACCTAATGAGAATTTT  
 CAGTCACTTTATAATAAGGAATTGCCTGTGCATATCTGTAATGTAATATCTCCTGAGAAGATTTATGTTT  
 AGTGGTTGTTAACTGAAAACCTACTTAATAGTTTAGAAGAAAAGATGATAGCTGCTTATGAAAACCTAAA  
 ATGGGAACCTGTTAAATGGGAAAATGATATGCACTGTGCTGTTAAGATCCAAGATAAAAAACAGTGGCGA  
 AGAGGCCAGATCATCAGAATGGTTACAGACACATTGGTAGAGGTCTTGCTGTATGATGTGGGTGTTGAAC  
 TAGTAGTGAATGTTGACTGTTTTAAGAAAACCTGAAGAAAATCTAAGACAATGGGAAGACTCTCTTTGGA  
 ATGTTCTCTGTTGACATAAGACCAGCTGGTGGAGTGACAAGTGAGCAACAGCTTGTGACTTGTGACTCTT  
 TCATTGTACCTGACTGGAGCTGTAGCAACTATAATCTTACAGGAAAACAACACAACATGGCCATTACCTG  
 TGAAAAATTTCTGCAGAGATGAAAAAGGAGAGCGTGTTGATGTTTCTAAATATTTGATTA AAAAGGGTTT  
 GGCTTTGAGAGAAAGGAGAATTAATAACTTAGATAACAGCCATTCTATCTGAGAAGTCTCTGGAAGTC  
 CCCCTGGAACAGGAAGATTCAAGTACTAAGTGTATTTAACTAACTTTGACCCTGACAAGAAAACCTG  
 CTGACATAATCAGTGAACAGAAAGTGTCTGAATTTCAAGGAGAAAATCTAGAACCAAGAACCACTAGAGG  
 GTATAAGCCACCAGCTATTCCTAACATGAACGTATTTGAGGCAACAGTCAGCTGTGTTGGTGATGATGGA  
 ACTATATTTGTAGTACCTAACTATCAGAAATTTGAGCTAATAAAAATGACAAATGAAATTCAAAGTAATT  
 TAAAAATGCCTTGGTCTTTTGGAGCCTTATTTCTGGA AAAAGGAGAAGCATGTGCAGTAAGAGGATCCGA  
 TACTCTGTGGTATCGTGGCAAGGTGATGGAGGTTGTAGGTGGCGCTGTGAGAGTACAATATTTAGATCAT  
 GGATTCAGTGA AAAGATTCCGCACTGCCATCTTACCCTATTTTGTCTGTATCCTGATATACCCAGTTTT  
 GTATTCCTTGTGACGTCCATAATACCACACCTGTTGGGAATGTCTGGCAACCAGATGCAATAGAAGTTCT  
 TCAACAACTGCTTTCAAAGAGACAGGTGGACATTACATTATGGAGTTACCTAAAAATCCATGGGAGAAA  
 TTGTCTATTCACCTCTATTTTGTGGAATGTCACTTTCTATTTTATGGCATACTATAAACTGTACTT  
 CTGAACATACTGAGGAGATGTTGAAAGAAAAACCAAGATCAGATCATGATAAAAAGTATGAAGAGGAACA  
 ATGGGAAAATAAGGTTTGAGGAATTGCTTTGCGCTGAAACAGACACTCCTCTTTTACCACCATATTTGTCT  
 TCATCTCTGCCTTCCCAGGAGAACTCTATGCTGTTCAAGTTAAGCACGTTGTCTCACCTAATGAAGTGT  
 ATATTTGCCTTGATTCTATAGAACTTCTAACCAGTCTAACCAGCATAGTGACACAGATGATAGTGGAGT  
 CAGCGGGGAATCAGAATCCGAGAGCCTTGATGAAGCACTGCAGAGGGTTAATAAGAAGGTAGAGGCGCTT  
 CCTCCTCTGACGGATTTTGAACAGAAATGCCTTGCTTGCAGAAATGATGATGGCTTATGGTATAGAG  
 CGAAGATTGTTGCCATTAAGAATTTAATCCTTTATCTATCTTAGTACAATTTGTTGATTATGGATCAAC  
 TGCAAAGCTGACATTAACAGACTGTGCCAAATTCCTTCTCATCTTATGCGGTATCCAGCTCGAGCCATA  
 AAGGTTCTCTTGGCAGGGTTTAACTCCCTTAAGGGATCTAGGGGAGACAAGAATACCATATTGCCCCA  
 AATGGAGCATGGAGGCACTGTGGCTATGATAGACTGTCTTCAAGGAAAACAACCTATGCTGTGTCCAT  
 GGCTCCAGCACCAGAACAGATAGTGACATTATATGACGATGAACAGCATCCAGTTTCATATGCCGTTGGTA  
 GAAATGGGGCTTGCAGATAAAGATGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


**Protein Sequence:** >RC230725 representing NM\_001184993  
Red=Cloning site Green=Tags(s)

MAAEASKTGPSRSSYQRMGRKSQPWGAEEIQCTRCRRVSRSSGHHCELQCGHAFCELCLLMTTEECTII  
CPDCEVATAVNTRQRYYPMAGYIKEDSIMKQLPKTIKNCSQDFKKTADQLTTGLERSASTDKTLLNSSA  
VMLDTNTAAEEIDEALNTAHSFEQLSIAGKALEHMQKQTIERERVIEVEKQFDQLLAFDSTRKKNLCE  
EFARTDDYLSNLIKAKSYIEEKNNLNAAMNIARALQLSPSLRTYCDLNQIIRTLQLTSDSELAQVSSP  
QLRNPRLSVNCSEIICMFNNMGKIEFRDSTKCYQENEIRQNVQKKYNNKELSCYDTPPLEKKKVD  
SVLTSEAPPPLQPETNDVHLEAKNFQPKQDVATASPKTIAVLPMGSSPDVIEEIIEDNVESSELVF  
VSHVIDPCHFYIRKYSQIKDAKLEKKVNEFCNRSSHLDPDILELGARIFVSSIKNGMWCRTITELIP  
IEGRNTRKPCSPTRLFVHEVALQIFMVDGENSEVLITGVVDTHVRPEHSAKQHIALNDLCLVLRKSEP  
YTEGLLDKIQLPAQPCSLKDIPQNSNEGWEAAKVEFLKMVNNKAVSMKVFREEDGVLIVDLQKPPNK  
ISSDMPVSLRDALVFMELAKFKSQSLRSHFEKNTTLHYHPPILPKEMTDVSVTVCHINSPGDFYQLIEG  
LDILFLLKTIIEEFYKSEGENLEILCPVQDQACVAKFEDGIWYRAKVIGLPGHQVEVYKVDGNTAKIT  
IKDVRKIKDEFNLNAPEKAICKLAYIEPYKRTMQWSKEAKEFEKAQDKFMTCSVIKILEDNVLLVELF  
DSLGAPEMTTTSINDQLVKEGLASYEIGYILKDNSQKHIEVWDPSPREEIISNEVHNLNPVSAKSLPNENF  
QSLYNKELPVHICNVISPEKIYVQWLLTENLLNSLEEKMIAAAYENSKWEPVKWENDMHCARKIQDNQWR  
RGQIIRMVDTLVEVLLYDVGVELVNVNDCRLKLEENLKTMGRLSLECSLVDIRPAGGSDKWTATACDCL  
SLYL TGAVATIIILQENNTTWPLPVKIFCRDEKGERVDVSKYL IKKGLALRERRINNLDNSHSLSEKSLE  
PLEQEDSVVNTNICTNFDPDKKTADIISEQKVSEFQEKILEPRTTTRYKPPAIPNMNVFEATVSCVGDDG  
TIFVVPKLSSEFELIKMTNEIQSNLKLGLLEPYFWKKGEACAVRGSDTLWYRGKVMVVGGAVRVQYLDH  
GFTEKIPQCHLYPILLYPDIPQFCIPQLHNTPVGNVWQPD AIEVLQQLSKRQVDIHIMELPKNPWEK  
LSIHL YFDGMSLSYFMAYYKYCTSEHTEMLKEKPRSDHDKKYEQQWEIRFEELL SAETDTPLLPPYLS  
SSLPSGELYAVQVKHVSPNEVYICLDSIETSNQSNQHSDDSGVSGESESESLDEALQRVNNKVEAL  
PPLTDFRTEMPPLAEYDDGLWYRAKIVAIKEFNPLSILVQFVDYGSTAKLTLNRLCQIPSHLMRYPARAI  
KVLLAGFKPPLRDLGETRIPYCPKWSMEALWAMIDCLQKGQLYAVSMAPAPEQIVTLYDDEQHPVHMPLV  
FMGLADKDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI


### Cloning Scheme:



ACCN: NM 001184993

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORF Size:</b>              | 4857 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_001184993.2</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq ORF:</b>            | 4860 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 56163                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Cytogenetics:</b>          | 13q12.12                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>MW:</b>                    | 184.7 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Gene Summary:</b>          | This gene is similar to a mouse gene that encodes a testis-specific protein containing a RING finger domain. Alternatively spliced transcript variants encoding different isoforms have been found. [provided by RefSeq, May 2010]                                                                                                                                                                                                                                                                                      |

## Product images:



Circular map for RC230725