

Product datasheet for **RC230725L4V**

RNF17 (NM_001184993) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RNF17 (NM_001184993) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF17
Synonyms:	Mmip-2; SPATA23; TDRD4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001184993
ORF Size:	4857 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230725).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001184993.1
RefSeq ORF:	4860 bp
Locus ID:	56163
Cytogenetics:	13q12.12
MW:	184.7 kDa
Gene Summary:	This gene is similar to a mouse gene that encodes a testis-specific protein containing a RING finger domain. Alternatively spliced transcript variants encoding different isoforms have been found. [provided by RefSeq, May 2010]


[View online »](#)