

Product datasheet for RC230777L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

glutathione S transferase Omega 1 (GSTO1) (NM_001191003) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: glutathione S transferase Omega 1 (GSTO1) (NM 001191003) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: glutathione S transferase Omega 1

Synonyms: GSTO 1-1; GSTTLp28; HEL-S-21; P28; SPG-R

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC230777).

Sequence:

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:






^{*} The last codon before the Stop codon of the ORF.





Plasmid Map:



ACCN: NM_001191003

ORF Size: 726 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001191003.1, NP 001177932.1</u>

RefSeq Size: 1070 bp RefSeq ORF: 642 bp



glutathione S transferase Omega 1 (GSTO1) (NM_001191003) Human Tagged Lenti ORF Clone – RC230777L4

 Locus ID:
 9446

 UniProt ID:
 P78417

 Cytogenetics:
 10q25.1

Protein Families: Druggable Genome

Protein Pathways: Drug metabolism - cytochrome P450, Glutathione metabolism, Metabolism of xenobiotics by

cytochrome P450

MW: 27.6 kDa

Gene Summary: The protein encoded by this gene is an omega class glutathione S-transferase (GST) with

glutathione-dependent thiol transferase and dehydroascorbate reductase activities. GSTs are involved in the metabolism of xenobiotics and carcinogens. The encoded protein acts as a homodimer and is found in the cytoplasm. Three transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Jul 2010]