

Product datasheet for RC230830L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MADH7 (SMAD7) (NM 001190821) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: MADH7 (SMAD7) (NM_001190821) Human Tagged ORF Clone Lentiviral Particle

Symbol: SMAD7

Synonyms: CRCS3; MADH7; MADH8

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_001190821

ORF Size: 1275 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC230830).

Sequence:

Cytogenetics:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This

clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 001190821.1, NP 001177750.1

18q21.1

 RefSeq Size:
 3085 bp

 RefSeq ORF:
 1278 bp

 Locus ID:
 4092

 UniProt ID:
 015105

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: TGF-beta signaling pathway





MW: 46.4 kDa

Gene Summary: The protein encoded by this gene is a nuclear protein that binds the E3 ubiquitin ligase

SMURF2. Upon binding, this complex translocates to the cytoplasm, where it interacts with TGF-beta receptor type-1 (TGFBR1), leading to the degradation of both the encoded protein and TGFBR1. Expression of this gene is induced by TGFBR1. Variations in this gene are a cause of susceptibility to colorectal cancer type 3 (CRCS3). Several transcript variants encoding

different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]