

Product datasheet for **RC230848L3V**

CHST15 (NM_014863) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CHST15 (NM_014863) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CHST15
Synonyms:	BRAG; GALNAC4S-6ST
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014863
ORF Size:	1518 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230848).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014863.3
RefSeq Size:	3580 bp
RefSeq ORF:	1521 bp
Locus ID:	51363
Cytogenetics:	10q26.13
Domains:	Sulfotransfer
Protein Families:	Transmembrane
Protein Pathways:	Chondroitin sulfate biosynthesis



[View online »](#)

MW: 58.1 kDa

Gene Summary: Chondroitin sulfate (CS) is a glycosaminoglycan which is an important structural component of the extracellular matrix and which links to proteins to form proteoglycans. Chondroitin sulfate E (CS-E) is an isomer of chondroitin sulfate in which the C-4 and C-6 hydroxyl groups are sulfated. This gene encodes a type II transmembrane glycoprotein that acts as a sulfotransferase to transfer sulfate to the C-6 hydroxyl group of chondroitin sulfate. This gene has also been identified as being co-expressed with RAG1 in B-cells and as potentially acting as a B-cell surface signaling receptor. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2012]