

Product datasheet for **RC230854**

ORC3L (ORC3) (NM_001197259) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORC3L (ORC3) (NM_001197259) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ORC3
Synonyms:	LAT; LATHEO; ORC3L
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC230854 representing NM_001197259
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTACGTCCTCGATGTCTAAGGGTTGCTTTGTTTTAAAGCCAACTCCAAAAAGAGAAAGATCTCTC
 TGCCAATAGAGGACTATTTAAACAAAGGAAAAATGAGCCTGAGGACAGTAAGCTTCGATTCGAACTTA
 TCAGTTGATATGGCAGCAGATGAAATCTGAAAATGAGCGACTACAAGAGGAATTAATAAAAACTTGTTT
 GACAATCTGATTGAATTTCTGCAAAAATCACATTCTGGATTCCAGAAGAATCAAGAGACTTGGCGGGT
 AAATAAACTCAGAGAAATCCAACTGCTGCTCTTGTCTTGGTGTGAATGTCACAGATCATGATTTGAC
 ATTCGGAAGTCTAACAGAGGCCCTTCAGAATAATGTCACACCATATGTAGTCTCATTGCAAGCTAAAGAT
 TGTCAGATATGAAACATTTTTGCAAAAGTTGATCTCACAGTTGATGGACTGCTGTGTAGATATAAAAT
 CCAAAGAGGAGGAAAGTGTTCACGTACCCAAAGAAAGACACATTATTCAATGGATTCACTTTCCAGTTG
 GTATATGACTGTCACACAGAAGACGGACCCAAAAATGCTAAGCAAAAAAGGACTACTTCTAGCCAAATGG
 CAGTCTCCTCCTGTTGTGCTTATCTTGAAGGATATGAAAAGCTTTGCCACAAAAGTACTACAAGACTTCA
 TAATTATCAGCAGTCAACATCTCCATGAATTTCCACTAATACTCATTTTTGGAATAGCCACATCTCCTAT
 TATCATCCACCGATTGCTTCTCATGCAGTATCATCTCTATTGTGCATAGAACTGTTCCAATCTTTGCT
 TGTAAGGAGCACCTGACTACGGTACTCGATAAGCTACTTCTTACAACCTCAGTTTCCCTTTAAAAATAATG
 AAAAGTATTACAGGTTCTGACCAACATCTTTTGTATCATGATTTCTCAGTTCAAACTTTATAAAAGG
 ACTTCAGTCTTCTCTATTAGAGCATTCTATTCCAGCCCTTAAGTGTCTGTGCTGTAATCTTCCAGAA
 GCCAAAAGAAGAATAAATTTTTATCAAATAATCAATGTGAAAACATCCGACGTCTACCATCTTTTAGGA
 GGTACGTGGAAAAGCAAGCTTCAGAAAAGCAAGTTGCGCTCTTGACCAATGAGAGATATTTGAAGGAGGA
 AACACAATTATTACTAGAAAACCTGCATGTTTATCATATGAATTACTTCTGGTTTTGAGATGTCTTCAT
 AAGTTCACCTCTTCTTCCCAAGTATCCACTAGGTCGACAGATCAGAGAGTTGACTGTACATGTTTAG
 AAAAGAACATATGGGATTCAGAGGAGTATGCATCAGTCTTGCAGCTGCTGAGGATGTTGGCAAAGGATGA
 ACTGATGACCATACTTGAGAAATGTTTCAAGGTTTTAAGTCTTATTGTGAAAACCACTTGGCAGCACA
 GCTAAGAGAATAGAGGAGTTCCTGGCCAGTTTCAGAGCCTCGATGAAACCAAAGAGGAAGAAGATGCTT
 CTGGGTACAGCCAAAGGGGCTTCAGAAGACAGACCTCTATCATCTTCAGAAGTCCTTATTGAAATGAA
 GGAGTTAAGAAGAAGTAAGAAGCAAACCAATTTGAAGTACTCAGAGAAAATGTTGTGAACCTCATTGAC
 TGTCTAGTGAGAGAATACCTTCTGCCTCCTGAGACACAGCCTCTCCATGAGGTGGTGTACTTCAGTGCTG
 CCCATGCCCTTCGTGAGCATTTAAATGCTGCTCCGCGAATTGCCCTCCATCTGCACTCAACAATCCTTA
 CTATTATCTCAAGAATGAAGCACTGAAAAGCGAAGAAGGCTGCATTCCGAATATCGCCCCAGACATCTGC
 ATAGCATACAACTGCACCTAGAGTGTAGCAGGCTCATCAACCTCGTGGACTGGTCAGAGGCTTTTGCAA
 CAGTTGTGACAGCTGCTGAAAAATGGATGCAAACTTCTGCAACCTCAGAAGAAAATGAATGAAATTTCCA
 TGCTCGGTTTATTAGAGCTGTTTCTGAAGTACTGAACTTTTAGGATTTATAAAACCTACCAACAGAAGACT
 GACCATGTGGCAAGACTAACATGGGGAGGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC230854 representing NM_001197259
 Red=Cloning site Green=Tags(s)

MATSSMSKGCVFVKPNSKKRKISLPIEDYFNKGKNEPEDSKLRFETYQLIWQQMKSENERLQEELNKNLFDNLI EFLQKSHSGFQKNSRDLGGQIKLREIPTAALVLGVNVDHDLTFGSLTEALQNNVTPYVVS LQAKDCPDMKHF LQKLISQLMDCCVDIKSKEEESVHVVTQRKTHYSMDLSSWYMTVTQKTDPKMLSKKRTTSSQWQSPPPVVILKDMESFATKVLQDFIIISSQHLHEFPLILIFGIATSPIIIHRLPHAVSSLLCIELFQSLSCKEHLTTVLDKLLLTTFPFKINEKVLQVLTNIFLYHDFSVQNFIKGLQLSLEHFYSQPLSVLCCNLPEAKRRINFLSNNQCENIRRLPSFRRYVEKQASEKQVALLTNERYLKEETQLLLENLHVYHMNYFLVLRCLHKFTSSLPKYPLGRQIRELYCTCLEKNIWDEEYASVLQLLRMLAKDELMTILEKCFKVFKSYCENHLGSTAKRIEEFLAQFQSLDEKKEEDASGSQPKGLQKTDLYHLQKSLEMKELRRSKKQTKFEVLRNVNFIDCLVREYLLPPETQPLHEVVYFSAHALREHLNAPRIALHTALNPPYYLKNEALKSEEGCIPNIAPDICIAYKLHLECSRLINLVDWSEAFATVVTAAEKMDANSATSEEMNEIIHARFIRAVSELELLGFIKPTKQKTDHVARL TWGGC


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001197259

ORF Size: 2136 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001197259.1](#), [NP_001184188.1](#)

RefSeq Size: 2488 bp

RefSeq ORF: 1707 bp

Locus ID: 23595

UniProt ID: [Q9UBD5](#)

Cytogenetics: 6q15

Protein Pathways: Cell cycle

MW: 82.3 kDa

Gene Summary: The origin recognition complex (ORC) is a highly conserved six subunits protein complex essential for the initiation of the DNA replication in eukaryotic cells. Studies in yeast demonstrated that ORC binds specifically to origins of replication and serves as a platform for the assembly of additional initiation factors such as Cdc6 and Mcm proteins. The protein encoded by this gene is a subunit of the ORC complex. Studies of a similar gene in *Drosophila* suggested a possible role of this protein in neuronal proliferation and olfactory memory. Alternatively spliced transcript variants encoding distinct isoforms have been reported for this gene. [provided by RefSeq, Jul 2008]