

Product datasheet for **RC230861**

Annexin VI (ANXA6) (NM_001193544) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Annexin VI (ANXA6) (NM_001193544) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Annexin VI
Synonyms:	ANX6; CBP68; CPB-II; p68; p70
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC230861 representing NM_001193544
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCCAAACCAGCACAGGGTGCCAAGTACCGGGGCTCCATCCATGACTTCCCAGGCTTTGACCCCAACC
 AGGATGCCGAGGCTCTGTACACTGCCATGAAGGGCTTTGGCAGTGACAAGGAGGCCATACTGGACATAAT
 CACCTCACGGAGCAACAGGCAGAGGCAGGAGGTCTGCCAGAGCTACAAGTCCCTCTACGGCAAGGACCTC
 ATTGCTGATTTAAAGTATGAATTGACGGCAAGTTGAACGGTTGATTGTGGCCTGATGAGGCCACCTG
 CCTATTGTGATGCCAAAGAAATTAAGATGCCATCTCGGGCATTGGCACTGATGAGAAGTGCCTCATTGA
 GATCTTGGCTTCCCGACCAATGAGCAGATGCACCAGCTGGTGGCAGCATACAAAGATGCCTACGAGCGG
 GACCTGGAGGCTGACATCATCGGCACACCTCTGGCCACTCCAGAAGATGCTTGTGGTCTGCTCCAGG
 GAACCAGGGAGGAGGATGACGTAGTGAGCGAGGACCTGGTACAACAGGATGTCCAGGACCTATACGAGGC
 AGGGAACTGAAATGGGGAACAGATGAAGCCCAGTTTACATCTTGGGAAATCGCAGCAAGCAGCAT
 CTTCCGTTGGTGTTCGATGAGTATCTGAAGACCACAGGGAAGCCGATTGAAGCCAGCATCCGAGGGGAGC
 TGCTGGGGACTTTGAGAAGCTAATGCTGGCCGTAGTGAAGTGTATCCGGAGCACCCCGGAATTTTTGC
 TGAAAGGCTCTTCAAGGCTATGAAGGGCTGGGGACTCGGGACAACCCCTGATCCGCATCATGGTCTCC
 CGTAGTGAGTTGGACATGCTCGACATTCGGGAGATCTTCCGGACCAAGTATGAGAAGTCCCTCTACAGCA
 TGATCAAGAAATGACACCTCTGGCGAGTACAAGAAGACTCTGCTGAAGCTGTCTGGGGGAGATGATGATGC
 TGCTGGCCAGTTCTTCCCGGAGGCAGCGCAGGTGGCCTATCAGATGTGGGAACTTAGTGCAGTGGCCGA
 GTAGAGCTGAAGGAACTGTGCGCCAGCCAATGACTTCAACCTGACGCAGATGCCAAAGCGCTGCGGA
 AAGCCATGAAGGGACTCGGGACTGACGAAGACACAATCATCGATATCATCAGCACCGCAGCAATGTCCA
 GCGGCAGCAGATCCGGCAGACCTTCAAGTCTCACTTTGGCCGGGACTTAATGACTGACCTGAAGTCTGAG
 ATCTCTGGAGACCTGGCAAGGCTGATTCTGGGGCTCATGATGCCACCGGCCATTACGATGCCAAGCAGT
 TGAAGAAGGCCATGGAGGGAGCCGGCACAGATGAAAAGGCTCTTATTGAAATCCTGGCCACTCGGACCAA
 TGCTGAAATCCGGGCCATCAATGAGGCCTATAAGGAGGACTATCACAAAGTCCCTGGAGGATGCTCTGAGC
 TCAGACACATCTGGCCACTTCAGGAGGATCCTCATTTCTCTGGCCACGGGGCATCGTGAGGAGGGAGGAG
 AAAACCTGGACCAGGCACGGGAAGATGCCAGGTGGCTGCTGAGATCTTGGAAATAGCAGACACACCTAG
 TGGAGACAAAACCTTCTTGGAGACACGTTTCATGACGATCCTGTGTACCCGGAGCTATCCGCACCTCCGG
 AGAGTCTTCCAGGAGTTCATCAAGATGACCAACTATGACGTGGAGCACACCATCAAGAAGGAGATGTCTG
 GGGATGTCAGGGATGCATTTGTGGCCATTGTTCAAAGTGTCAAGAACAAGCCTCTCTTTTGGCCGACAA
 ACTTTACAAATCCATGAAGGGTCTGGCACAGATGAGAAGACTCTGACCAGGATCATGGTATCCCGCAGT
 GAGATTGACCTGCTCAACATCCGGAGGGAATTCATTGAGAAATATGACAAGTCTCTCCACCAAGCCATTG
 AGGGTGACACCTCCGGAGACTTCTGAAGGCCTTGTGGCTCTCTGTGGTGGTGAGGAC

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC230861 representing NM_001193544
Red=Cloning site Green=Tags(s)

```
MAKPAQGAKYRGSIHDFPGFDPNQDAEALYTAMKGFSGDKEAILDITSRNRQRQEVQSYKSLYGKDL
IADLKYELTGKFERLIVGLMRPPAYCDAKEIKDAISGIGTDEKCLIEILASRTNEQMHQLVAAYKDAYER
DLEADIIIGDTSGHFQKMLVLLQGTREEDDVVSEDLVQQDVQDLYEAGELKWGTDEAQFIYILGNRSKQH
LRLVDFEYLKTTGKPIEASIRGELSGDFEKLMLAVVKCIRSTPEYFAERLFAKMKGLGTRDNTLIRIMVS
RSELDMLDIREIFRTKYEKSLYSMIKNDTSGEYKTLKLSGGDDDAAGQFFPEAAQVAYQMWELSAVAR
VELKGTVRPANDFNPDADAKALRKAMKGLGTDEDIIDIIITHRSNVQRQQIRQTFKSHFGRDMLTKLSE
ISGDLARLILGLMPPAHYDAKQLKKAMEGAGTDEKALIEILARTNAEIRAINAEYKEDYHKSLEDALS
SDTSGHFRRILISLATGHREEGGENLDQAREDAQVAEILEIADTPSGDKTSLETRFMTILCTRSYPHLR
RVFQEFIKMTNYDVEHTIKKEMSGDVRDAFVAIVQSVKNKPLFFADKLYKSMKGAGTDEKTLTRIMVSR
EIDLLNIRREFIEKYDKSLHQAIEGDTSGDFLKALLALCGGED
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001193544

ORF Size: 2019 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001193544.1](#), [NP_001180473.1](#)

RefSeq Size: 2879 bp

RefSeq ORF: 1926 bp

Locus ID: 309


UniProt ID: [P08133](#)

Cytogenetics: 5q33.1

MW: 75.9 kDa

Gene Summary: Annexin VI belongs to a family of calcium-dependent membrane and phospholipid binding proteins. Several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. The annexin VI gene is approximately 60 kbp long and contains 26 exons. It encodes a protein of about 68 kDa that consists of eight 68-amino acid repeats separated by linking sequences of variable lengths. It is highly similar to human annexins I and II sequences, each of which contain four such repeats. Annexin VI has been implicated in mediating the endosome aggregation and vesicle fusion in secreting epithelia during exocytosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010]

Product images:



Circular map for RC230861

