

Product datasheet for **RC230932A1V**

Human **FAM240A (NM_001195442) AAV Particle**

Product data:

Product Type:	AAV Particles
Product Name:	Human FAM240A (NM_001195442) AAV Particle
Tag:	Myc-DDK
Symbol:	FAM240A
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC230932 representing NM_001195442 Red =Cloning site Blue =ORF Green =Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC ATGCTTCACATCCCTTCATGCGGTTGTTCTCTGGGATGAACAATCAATACACCCGTCGGGAGGTCTTCTGCCGGAACACCTGCCATGATCTCAAGCATTCTGGGAAAGGGAAATTGGCAAACAGACTTACTACCGAGAATCAGAGGAACGTCGTCGGAAGAAGCGCACTGAGAAAGCTCAGAGAAGAATGGAAGCAGAGGCTGGAAACAAGCTGAGGCTCGGGAACAATCCAGAGGACACTGAAAAGAGGACAAATGTTGGC ACGCGT ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC230932 representing NM_001195442 Red =Cloning site Green =Tags(s) MSSHPPMRLFSGMNNQYTRREVF CRNTCHDLKHFWEREIGKQTYRESEERRLGRSALRKLREEWKQRLETKLRLRNNPEDTEKRTNVG TR TRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001195442
ORF Size:	267 bp
Buffer:	PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_001195442.1](#), [NP_001182371.1](#)

RefSeq ORF: 252 bp

Locus ID: 100132146

UniProt ID: [A0A1B0GVK7](#)

Cytogenetics: 3p21.31

MW: 11.5 kDa