

## Product datasheet for **RC230963L4V**

### AF10 (MLLT10) (NM\_001195627) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	AF10 (MLLT10) (NM_001195627) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MLLT10
Synonyms:	AF10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001195627
ORF Size:	378 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230963).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001195627.1</a> , <a href="#">NP_001182556.1</a>
RefSeq ORF:	381 bp
Locus ID:	8028
UniProt ID:	<a href="#">P55197</a>
Cytogenetics:	10p12.31
Protein Families:	Druggable Genome, Transcription Factors
MW:	14.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a transcription factor and has been identified as a partner gene involved in several chromosomal rearrangements resulting in various leukemias. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2010]