

## Product datasheet for **RC230984L2V**

### Gremlin 1 (GREM1) (NM\_001191323) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Gremlin 1 (GREM1) (NM_001191323) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Gremlin 1  |
| Synonyms:                 | C15DUPq; CKTSF1B1; CRAC1; CRCS4; DAND2; DRM; DUP15q; GREMLIN; HMPS; HMPS1; IHG-2; MPSH; PIG2   |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-mGFP (PS100071)   |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001191323   |
| ORF Size:                 | 429 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC230984).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001191323.1</a>   |
| RefSeq Size:              | 4015 bp  |
| RefSeq ORF:               | 432 bp   |
| Locus ID:                 | 26585  |
| UniProt ID:               | <a href="#">O60565</a>   |
| Cytogenetics:             | 15q13.3  |
| Protein Families:         | ES Cell Differentiation/IPS, Secreted Protein  |



[View online »](#)

**MW:** 16.3 kDa

**Gene Summary:** This gene encodes a member of the BMP (bone morphogenic protein) antagonist family. Like BMPs, BMP antagonists contain cystine knots and typically form homo- and heterodimers. The CAN (cerberus and dan) subfamily of BMP antagonists, to which this gene belongs, is characterized by a C-terminal cystine knot with an eight-membered ring. The antagonistic effect of the secreted glycosylated protein encoded by this gene is likely due to its direct binding to BMP proteins. As an antagonist of BMP, this gene may play a role in regulating organogenesis, body patterning, and tissue differentiation. In mouse, this protein has been shown to relay the sonic hedgehog (SHH) signal from the polarizing region to the apical ectodermal ridge during limb bud outgrowth. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2010]