

## Product datasheet for RC230986A1V

### Human KIAA1522 (NM\_001198973) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human KIAA1522 (NM\_001198973) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** KIAA1522  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC230986 representing NM\_001198973  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGTGGTGTTTCGTTGGCCGCCCTCCCGGCGCTCCTAGGGCTGTTTAAGAAGAAGGGCTCTGCCAAGG  
 CTGAGAATGACAAACATCTAAGTGTAGGGCCTGGCCAGGGGCCAGGGTCTGCAGTGGATGAGCACCAGGA  
 CAACGTCTTCTTCCAGTGGGCGACCCCCACCTGGAAGAGCTGCACACTCAGGCCAGGAGGGGCTC  
 CGCTCCCTACAACACCAAGAGAAACAGAACTGAACAAGGGTGGCTGGGACCATGGAGACACCCAGAGTA  
 TCCAGTCTCCCGGACGGGGCCGGATGAAGACAACATCTCCTTCTGCAGTCAGACCACATCCTACGTGGC  
 TGAGAGCTCCACAGCAGAGGACGCGCTCCATCCGCTCGGAGATGATCCAGCGCAAAGACTCACAGAAA  
 GAGCTGGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC230986 representing NM\_001198973  
 Red=Cloning site Green=Tags(s)

MVVFVGRRLPALLGLFKKGSKAENDKHLVGPQGPQSAVDEHQDNVFFPSGRPPHLEELHTQAQEGL  
 RSLQHQEKQKLNKGGWDHGDTSIQSSRTGPEEDNISFCSQTSYVAESSTAEDALSIRSEMIQRKDSQK  
 ELA

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001198973



[View online »](#)

<b>ORF Size:</b>	429 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001198973.1</a></u>
<b>RefSeq ORF:</b>	432 bp
<b>Locus ID:</b>	57648
<b>UniProt ID:</b>	<u><a href="#">Q9P206</a></u>
<b>Cytogenetics:</b>	1p35.1
<b>MW:</b>	16.1 kDa