

## Product datasheet for **RC230997A1V**

### Human CLEC2D (NM\_001197317) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human CLEC2D (NM\_001197317) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** CLEC2D  
**Synonyms:** CLAX; LLT1; OCIL  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC230997 representing NM\_001197317  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCATGACAGTAACAATGTGGAGAAAGACATTACACCATCTGAATTGCCTGCAAACCCAGCAATAAGAG  
CTAACTGCCATCAAGAGCCATCAGTATGTCTTCAAGCTGCATGCCAGAAAGCTGGATTGGTTTTCAAAG  
AAAGTGTCTATTTTTCTGATGACACCAAGAAGCTGGACATCAAGTCAGAGGTTTTGTGACTCACAAGAT  
GCTGATCTTGCTCAGGTTGAAAGCTTCCAGGAAGTGAATTTCTGTTGAGATATAAAGGCCCATCTGATC  
ACTGGATTGGGCTGAGCAGAGAACAAGCCAACCATGGAATGGATAAAATGGTACTGAATGGACAAGACA  
GTTTTCTATCCTGGGAGCAGGAGAGTGTGCCTATTTGAATGACAAAGGTGCCAGTAGTCCAGGCACTAC  
ACAGAGAGGAAGTGGATTTGTTCCAATCAGATATACATGTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC230997 representing NM\_001197317  
Red=Cloning site Green=Tags(s)

MHDSNNVEKDITPSEL PANPAIRANCHQEPSVCLQAACPESWIGFQRKCFYFSDDTKNWTSSQRFCDSDQ  
ADLAQVESFQELNLLRYKGPSDHWIGLSREQQPWKINGTEWTRQFPILGAGECAYLNDKGASSARHY  
TERKWICKSIDIHV

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2



[View online »](#)

<b>ACCN:</b>	NM_001197317
<b>ORF Size:</b>	462 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001197317.1</a></u>
<b>RefSeq ORF:</b>	465 bp
<b>Locus ID:</b>	29121
<b>UniProt ID:</b>	<u><a href="#">Q9UHP7</a></u>
<b>Cytogenetics:</b>	12p13.31
<b>MW:</b>	18.2 kDa