

# **Product datasheet for RC231009**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## p16INK4A (CDKN2A) (NM 001195132) Human Tagged ORF Clone

#### **Product data:**

**Product Type: Expression Plasmids** 

**Product Name:** p16INK4A (CDKN2A) (NM\_001195132) Human Tagged ORF Clone

Tag: Myc-DDK Symbol: p16INK4A

Synonyms: ARF; CDK4I; CDKN2; CMM2; INK4; INK4A; MLM; MTS-1; MTS1; P14; P14ARF; P16; P16-INK4A;

P16INK4; P16INK4A; P19; P19ARF; TP16

Mammalian Cell

Selection:

Sequence:

Neomycin

Vector: pCMV6-Entry (PS100001) E. coli Selection: Kanamycin (25 ug/mL)

**ORF Nucleotide** 

>RC231009 representing NM\_001195132 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCGCGATCGCC

ATGGAGCCGGCGGGGAGCAGCATGGAGCCTTCGGCTGACTGGCCACGGCCGCCGGGCCCGGGGTC GGGTAGAGGAGGTGCGGGCGCTGCTGGAGGCGGGGGGCGCTGCCCAACGCACCGAATAGTTACGGTCGGAG GCCGATCCAGGTCATGATGATGGGCAGCGCCCGAGTGGCGGAGCTGCTGCTGCTCCACGGCGCGGAGCCC AACTGCGCCGACCCCGCCACTCTCACCCGACCCGTGCACGACGCTGCCCGGGAGGGCTTCCTGGACACGC TGGTGGTGCTGCACCGGGCCGGGCGCGGCTGGACGTGCCCGATGCCTGGGGCCGTCTGCCCGTGGACCT AACCATGCCCGCATAGATGCCGCGGAAGGTCCCTCAGAAATGATCGGAAACCATTTGTGGGTTTGTAGAA

**GCAGGCATGCG** 

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT

ACAAGGATGACGACGATAAGGTTTAA

>RC231009 representing NM\_001195132 **Protein Sequence:** 

Red=Cloning site Green=Tags(s)

MEPAAGSSMEPSADWLATAAARGRVEEVRALLEAGALPNAPNSYGRRPIQVMMMGSARVAELLLLHGAEP NCADPATLTRPVHDAAREGFLDTLVVLHRAGARLDVRDAWGRLPVDLAEELGHRDVARYLRAAAGGTRGS

NHARIDAAEGPSEMIGNHLWVCRSRHA

**TRTRPL**EQKLISEEDLAANDILDYKDDDDK**V** 



#### p16INK4A (CDKN2A) (NM\_001195132) Human Tagged ORF Clone - RC231009


https://cdn.origene.com/chromatograms/ja1865 a10.zip **Chromatograms:** 

**Restriction Sites:** 

Sgfl-Mlul

**Cloning Scheme:** 





ACCN: NM\_001195132

**ORF Size:** 501 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 001195132.2

RefSeq ORF: 504 bp Locus ID: 1029 UniProt ID: P42771



**Cytogenetics:** 9p21.3

**Protein Families:** Druggable Genome

Protein Pathways: Bladder cancer, Cell cycle, Chronic myeloid leukemia, Glioma, Melanoma, Non-small cell lung


cancer, p53 signaling pathway, Pancreatic cancer, Pathways in cancer

MW: 18.3 kDa

**Gene Summary:** This gene generates several transcript variants which differ in their first exons. At least three


alternatively spliced variants encoding distinct proteins have been reported, two of which encode structurally related isoforms known to function as inhibitors of CDK4 kinase. The remaining transcript includes an alternate first exon located 20 Kb upstream of the remainder of the gene; this transcript contains an alternate open reading frame (ARF) that specifies a protein which is structurally unrelated to the products of the other variants. This ARF product functions as a stabilizer of the tumor suppressor protein p53 as it can interact with, and sequester, the E3 ubiquitin-protein ligase MDM2, a protein responsible for the degradation of p53. In spite of the structural and functional differences, the CDK inhibitor isoforms and the ARF product encoded by this gene, through the regulatory roles of CDK4 and p53 in cell cycle G1 progression, share a common functionality in cell cycle G1 control. This gene is frequently mutated or deleted in a wide variety of tumors, and is known to be an important tumor suppressor gene. [provided by RefSeq, Sep 2012]

### **Product images:**



Circular map for RC231009





Western blot validation of overexpression lysate (Cat# [LY434008]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC231009 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).