

Product datasheet for **RC231010**

C16orf5 (CDIP1) (NM_001199055) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C16orf5 (CDIP1) (NM_001199055) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: CDIP1
Synonyms: C16orf5; CDIP; I1; LITAFL
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC231010 representing NM_001199055
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

**ATGTCCAGCGAGCCTCCCCTCCTTATCCTGGGGCCCCACAGCCCCACTTCTGGAAGAGAAAAGTGGAG
CCCCGCCACCCAGGCCGTTCCCTCCAGCTGTGATGCAGCCCCCTCCAGGCATGCCACTGCCCTCTG
GGACATTGGCCCCCACCCTATGAGCCCGGGTCACCCAATGCCCCAGCCTGGCTTCATCCACCACAC
ATGAGTGCAGATGGCACCTACATGCCTCCGGGAGCTGCCACCACGGTGACAGTGTGCAGGGAGAGATCT
TTGAGGGAGCGCCTGTGCAGACGGTGTGTCCCACTGCCAGCAGCCATCACCACCAAGATCTCCTACGA
GATTGGCTTGATGAATTCGTGCTGGGTTTCTTCTGTTGCTTCATGGGATGTGATCTGGGCTGTGCCTG
ATCCCCTGCCTCATCAATGACTTCAAGGATGTGACGCACACATGCCCCAGCTGCAAAGCCTACATCTACA
CGTACAAGCGCCTGTGC**

**ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**

Protein Sequence: >RC231010 representing NM_001199055
Red=Cloning site Green=Tags(s)


MSSEPPPPYPGGPTAPLLEEKSGAPPTGRSSPAVMQPPPGMPLPPADIGPPPYEPPGHMPQPQGFIPPH
MSADGTYMPPGAATTVTLQGEIFEGAPVQTVCPHCQQAITTKISYEIGLMNFVLGFFCCFMGCDLGCCL
IPCLINDFKDVTHTCPSCKAYIYTKRLC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:

Plasmid Map:


ACCN: NM_001199055

ORF Size: 507 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199055.2</u>
RefSeq ORF:	510 bp
Locus ID:	29965
UniProt ID:	<u>Q9H305</u>
Cytogenetics:	16p13.3
Protein Families:	Transmembrane
MW:	18.4 kDa
Gene Summary:	Acts as an important p53/TP53-apoptotic effector. Regulates TNF-alpha-mediated apoptosis in a p53/TP53-dependent manner.[UniProtKB/Swiss-Prot Function]