

## Product datasheet for **RC231012L4V**

### DHRS4L2 (NM\_001193635) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DHRS4L2 (NM_001193635) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DHRS4L2
Synonyms:	SDR25C3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001193635
ORF Size:	510 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231012).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001193635.1</a> , <a href="#">NP_001180564.1</a>
RefSeq ORF:	513 bp
Locus ID:	317749
Cytogenetics:	14q11.2
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Retinol metabolism
MW:	18.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the short chain dehydrogenase reductase family. The encoded protein may be an NADPH dependent retinol oxidoreductase. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2010]