

Product datasheet for RC231017A1V

Human TMEM238 (NM_001190764) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TMEM238 (NM_001190764) AAV Particle
Tag: Myc-DDK
Symbol: TMEM238
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC231017 representing NM_001190764
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGCGGCCAGCGGTGTGCGCCTCGCAGGGGAGCCCGCGGGTGACCGTCCGCGCCGGCGGCCG
 CGCCAGCACCCGCGCGGCCCTGGGCCGCTGCCGGATGGCGTGCTGCTGGCCGTGGCGCTGGATGTGGC
 GGGCATGGCGCGCTGCTGACCGCGTGTTCGCGCAGCTGCAGGTGCGCGGCCGCGACTTCGGGGACCTG
 CTCATCTACTCGGCGCGCTGCTGGTGTCTGAGCCTGCTGGGCTGGATCCTCTGGTACACCGCAACA
 TCGAGATCTCGCGCCAGGAGCTGGAGCGCGACTACGGCCTGCGGCCCTCGGCGCTCGCCGCCTGGCGCG
 CAAGCTCTCCCGCGCTGGTGGCGCCCCGCCCGCGGGCCAGCGCCCCGCGCCCGCTCCCGGAGAGCG
 CGCCGAGCCGCCCGCGCCCCCGCGCCCCGCCCGGCTCCCGCGCGTGCAGCTCGCCACGC
 TCGAGGCCGGGCCCGGGCGGCGGGCGCGGCAGCGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231017 representing NM_001190764
 Red=Cloning site Green=Tags(s)
 MAAAPAVCASQGSPPGAPSAPAAAPAPAAGLGRCRMALLLAVALDVAGMAALLTGVFQQLQVRGRDFGDL
 LIYSGALLVFLSLLGWILWYTGNIETSRQELERDYGLRPSALARLARKLSRRWSAPAAAGQRPAPGSRRA
 RRAARAPPPPAAGSRRVRLQLATLEAGPGAAGAGSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001190764



[View online »](#)

ORF Size:	528 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001190764.1
RefSeq ORF:	531 bp
Locus ID:	388564
UniProt ID:	C9JI98
Cytogenetics:	19q13.42
MW:	18.5 kDa