

Product datasheet for RC231050A1V

Human KLRF2 (NM_001190765) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: KLRF2

Synonyms: NKp65

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC231050 representing NM_001190765
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGAATGAAGATGGGTATATGACGCTGAGTTTCAAGAATCGTTGTAATCGAAGCAGAAATCTAAAG
ATTTCTCCCTATATCCACAATATTATTGTCTTCTGCTCATATTTGGATGCATTGTGATCCTTATATTCAT
TATGACAGGGATTGACCTGAAGTTCTGGCATAAAAAATGGATTTCTCCAGAATGTAACGTCAGCAGT
CTATCAGGACACAATTACTTGTGCCAAATGACTGGCTGTTGAACGAAGGGAATGTTACTGGTTTTCAA
CTTCTTTTAAACGTGGAAGAGAGTCAACGTGATTGTACACAGCTACAGGCACATTTACTGGTGATTCA
AAATTTGGATGAGCTGGAGTTCATACAGAACAGTTTAAACCTGGACATTTTGGTTGGATTGGACTATAT
GTTACATTCCAAGGGAACCTATGGATGTGGATAGATGAACACTTTTTAGTTCCAGAATTGTTTTCAAGTGA
TTGGACCAACTGATGACAGGAGCTGTGCCGTTATCACAGGAACTGGGTGATTCTGAAGACTGTAGCTC
CACATTTAAGGGCATTGCCAGAGAGATGCGATCTTGACGCACAATGGAACAGTGGTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231050 representing NM_001190765
Red=Cloning site Green=Tags(s)

MENEDGYMTLSFKNRCKSKQSKDFSLYPQYYCLLLIFGCIVILIFIMTGIDLKFWHKKMDFSQNVNVSS
LSGHNYLCPNDWLLNEGKCYWFSTSFKTWKESQRDCTQLQAHLLVIQNLDELEFIQNSLKPGHFGWIGLY
VTFQGNLWMWIDEHFLVPELFSVIGPTDDRSCAVITGNWVYSEDCSSTFKGICQRDAILTHNGTSGV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_001190765
ORF Size:	621 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001190765.1 , NP_001177694.1
RefSeq ORF:	624 bp
Locus ID:	100431172
UniProt ID:	D3W0D1
Cytogenetics:	12p13.31
MW:	24.5 kDa