

Product datasheet for **RC231050L3V**

KLRF2 (NM_001190765) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	KLRF2
Synonyms:	NKp65
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001190765
ORF Size:	621 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC231050).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001190765.1](#), [NP_001177694.1](#)

RefSeq ORF: 624 bp

Locus ID: 100431172

UniProt ID: [D3W0D1](#)

Cytogenetics: 12p13.31

MW: 24.5 kDa



Gene Summary:

C-type lectin-like receptor involved in natural killer cell mediated cytotoxicity and cytokine secretion in keratinocytes via its interaction with CLEC2A.[UniProtKB/Swiss-Prot Function]