

## Product datasheet for RC231057L4

### DKKL1 (NM\_001197301) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DKKL1 (NM_001197301) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	DKKL1
Synonyms:	CT34; SGY; SGY-1; SGY1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231057).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



UniProt ID: [Q9UK85](#)

Cytogenetics: 19q13.33

Protein Families: ES Cell Differentiation/IPS, Secreted Protein, Transmembrane

MW: 24.2 kDa

**Gene Summary:** The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Oct 2010]