

## Product datasheet for **RC231057L4V**

### **DKKL1 (NM\_001197301) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	DKKL1 (NM_001197301) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DKKL1
Synonyms:	CT34; SGY; SGY-1; SGY1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001197301
ORF Size:	633 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231057).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001197301.1</a>
RefSeq ORF:	636 bp
Locus ID:	27120
UniProt ID:	<a href="#">Q9UK85</a>
Cytogenetics:	19q13.33
Protein Families:	ES Cell Differentiation/IPS, Secreted Protein, Transmembrane
MW:	24.2 kDa



[View online »](#)

**Gene Summary:**

The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Oct 2010]