

Product datasheet for **RC231065L4V**

Apoptosis repressor with CARD (NOL3) (NM_001185057) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apoptosis repressor with CARD (NOL3) (NM_001185057) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Apoptosis repressor with CARD
Synonyms:	ARC; FCM; MYOCL1; MYP; NOP; NOP30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001185057
ORF Size:	657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231065).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001185057.1
RefSeq ORF:	660 bp
Locus ID:	8996
UniProt ID:	O60936
Cytogenetics:	16q22.1
MW:	24.8 kDa



[View online »](#)

Gene Summary:

This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]