

Product datasheet for **RC231110L4V**

EEF1D (NM_001195203) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | EEF1D (NM_001195203) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | EEF1D |
| Synonyms: | EF-1D; EF1D; FP1047 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001195203 |
| ORF Size: | 786 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC231110). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001195203.1 |
| RefSeq ORF: | 789 bp |
| Locus ID: | 1936 |
| UniProt ID: | P29692 |
| Cytogenetics: | 8q24.3 |
| MW: | 29.5 kDa |



[View online »](#)

Gene Summary:

This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit, delta, functions as guanine nucleotide exchange factor. It is reported that following HIV-1 infection, this subunit interacts with HIV-1 Tat. This interaction results in repression of translation of host cell proteins and enhanced translation of viral proteins. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene. Related pseudogenes have been defined on chromosomes 1, 6, 7, 9, 11, 13, 17, 19.[provided by RefSeq, Aug 2010]