

## Product datasheet for **RC231124L4V**

### **RNF138 (NM\_001191324) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

|                                  |                                   |
|----------------------------------|-----------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles              |
| <b>Symbol:</b>                   | RNF138                            |
| <b>Synonyms:</b>                 | hNARF; HSD-4; NARF; STRIN         |
| <b>Mammalian Cell Selection:</b> | Puromycin                         |
| <b>Vector:</b>                   | pLenti-C-mGFP-P2A-Puro (PS100093) |
| <b>Tag:</b>                      | mGFP                              |
| <b>ACCN:</b>                     | NM_001191324                      |
| <b>ORF Size:</b>                 | 846 bp                            |

**ORF Nucleotide Sequence:** The ORF insert of this clone is exactly the same as(RC231124).

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                      |                                |
|----------------------|--------------------------------|
| <b>RefSeq:</b>       | <a href="#">NM_001191324.1</a> |
| <b>RefSeq ORF:</b>   | 738 bp                         |
| <b>Locus ID:</b>     | 51444                          |
| <b>UniProt ID:</b>   | <a href="#">Q8WVD3</a>         |
| <b>Cytogenetics:</b> | 18q12.1                        |
| <b>MW:</b>           | 32.36 kDa                      |



**Gene Summary:**

The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]