

Product datasheet for **RC231128A1V**

Human **LRRC72 (NM_001195280)** AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LRRC72 (NM_001195280) AAV Particle
Tag:	Myc-DDK
Symbol:	LRRC72
Synonyms:	CFAP279; FAP279
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC231128 representing NM_001195280 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCCTGGGACCCGAACCCCGTGCCCCGTACCTTGCGATGCTGGCGCCTACGGAGGGCATCCGAACTG
CCCTACAGAGCAGTCGCCGGGCAGTTGAAGATCAGCTAAAGATATGTGGCCACAGGAGGGATGCTGATGT
CTTTGAGCTGTTCTTTCTAAAAAGGAAGTACAGAGGTCATTGATCTTTCTAGGTTTAAAAAATTTAAAA
TACTTATGGCTTCATCATAACAAGCTCCATGGAATAACATTTCTAACTAGAACTATTGTCTGACAGAAC
TATATCTAAACAACAATGCAATATTTGAGATAGAAGGTCTGCATTATTTGCCTTCATTGCATATACTCCT
GCTACACCACAATGAGCTAACCAACATTGATGCAACAGTGAAGGAATTAAGGGAATGCTAAATCTGAAG
ATCCTAAGTCTATACCAAAATCCTTTGTGCAATATAACCTGTATCGTTTATATATCATCTACCACTTC
CAGGAGTGGAGCTGCTTGACCGAAATCAAGTTACAGAAAAAGAAAGATCAATGATTACCATTTTAA
CCATAAAAGGCTCATATCGTTCAATCAATAGCATTGCGAGGAAAAGTGGATGCTTCATGGGATCCTAA
TCACCATTTAAGCAAAACCAGCCAGAGAGTACCTTCAGATTTTGCATTGCAATAATGTAGACAAAA
CTGTGCTTGATGACCCAGAAGATGCTGTTTTGTGAGGTCCATGAAGAGATCAGTGATGACTTTGACCTC
TATGAAGTGGGACACAGTTCCAACACGAGAGGAAAGGTACCTTGAAGAGGAAGGCACAGAAACAGCTCAA
ATGCTCACAGTTACACTGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC231128 representing NM_001195280 Red=Cloning site Green=Tags(s) MSWDPNPVPTLRCLRWLRRASETALQSSRAVEDQLKICGHRRDADVFEFLSKKELTEVIDLSRFKKLK YLWLHHNKLHGITFLTRNYCLTELYLNNAIFEIEGLHYLPSLHILLHHNELTNIDATVKELKGMNLK ILSLYQNPLCQYNLYRLYIIYHLPGVELLDRNQVTEKERRSMITIFNHKKAHIVQSIAFGGKVDASWDPK SPFKQKPAQRVPSDFAFANNVDKTVLDDPEDAVFVRSMKRSVMTLTSMNWDTPPTREERYLEEEGTETAQ MLTVTLR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001195280
ORF Size:	861 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001195280.1</u>
RefSeq ORF:	864 bp
Locus ID:	100506049
UniProt ID:	<u>A6NJI9</u>
Cytogenetics:	7p21.1
MW:	34.1 kDa