

## Product datasheet for **RC231133L3V**

### **AKR1D1 (NM\_001190907) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	AKR1D1 (NM_001190907) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AKR1D1
Synonyms:	3o5bred; CBAS2; SRD5B1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001190907
ORF Size:	870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231133).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001190907.1</a>
RefSeq ORF:	873 bp
Locus ID:	6718
UniProt ID:	<a href="#">P51857</a>
Cytogenetics:	7q33
Protein Families:	Druggable Genome
Protein Pathways:	Androgen and estrogen metabolism, C21-Steroid hormone metabolism, Metabolic pathways, Primary bile acid biosynthesis



[View online »](#)

**MW:** 33.2 kDa

**Gene Summary:** The enzyme encoded by this gene is responsible for the catalysis of the 5-beta-reduction of bile acid intermediates and steroid hormones carrying a delta(4)-3-one structure. Deficiency of this enzyme may contribute to hepatic dysfunction. Three transcript variants encoding different isoforms have been found for this gene. Other variants may be present, but their full-length natures have not been determined yet. [provided by RefSeq, Jul 2010]