

Product datasheet for **RC231158A1V**

Human DAZL (NM_001190811) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DAZL (NM_001190811) AAV Particle
Tag:	Myc-DDK
Symbol:	DAZL
Synonyms:	DAZH; DAZL1; DAZLA; SPGYLA
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC231158 representing NM_001190811 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCTCCCTCGTGTGGCGGCGACAGAAAAGCTCGCCTGACGCCATCTTTGCCGCACGAGTCTACTG
CAAATCCTGAACTCCAACTCAACCATCTCCAGAGAGGCCAGCACCCAGTCTCATCAGCTGCAACCAG
CCAAGGCTATATTTACCAGAAGGCAAAATCATGCCAAACACTGTTTTGTTGGAGGAATTGATGTTAGG
ATGGATGAACTGAGATTAGAAGCTTCTTTGCTAGATATGGTTCAGTGAAAGAAGTGAAGATAATCACTG
ATCGAACTGGTGTGTCCAAAGGCTATGGATTTGTTTCATTTTTTAATGACGTGGATGTCAGAAGATAGT
AGAATCAGAGATAAATTTCCATGGTAAAAAGCTGAAGCTGGGCCCTGCAATCAGGAAACAAAATTTATGT
GCTTATCATGTGCAGCCACGTCTTTGGTTTTTAATCATCCTCCTCCACCACAGTTTCAGAATGTCTGGA
CTAATCCAAACACTGAACTTATATGCAGCCCAACACGATGAATCCTATAACTCAGTATGTTCAAGC
ATATCCTACTTACCCAAATTCACCAGTTCAGGTCATCACTGGATATCAGTTGCCTGTATATAATTATCAG
ATGCCACCACAGTGGCTGTTGGGGAGCAAAGGAGCTATGTTGTACCTCCGGCTTATTCAGCTGTAACT
ACCACTGTAATGAAGTTGATCCAGGAGCTGAAGTTGTGCCAAATGAATGTTCAAGTTCATGAAGCTACTCC
ACCCTCTGGAATGGCCACAAAAGAAATCTGTGGACCGAAGCATACAAACGGTGGTATCTTGTCTGTTT
AATCCAGAGAACAGACTGAGAACTCTGTTGTTACTCAAGATGACTACTTCAAGGATAAAAGAGTGCATC
ACTTTAGAAGAAGTCGGGCAATGCTTAAATCTGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC231158 representing NM_001190811 Red=Cloning site Green=Tags(s) MAAPSCGGDRKARLTPSLPHESTANPETPNSTISREASTQSSSAATSQGYILPEGKIMPNTVFVGGIDVR MDETEIRSFFARYGSVKEVKIITDRTGVSKGYGFVSFFNDVDVQKIVESQINFHGKKLGLPAIRKQNL AYHVQPRPLVFNHPPPPQFQNVWTNPNTETYMQPTTTMNPITQYVQAYPTYPNSPVQVITGYQLPVYNYQ MPPQWPVGEQRSYVPPAYSAVNYHCNEVDPGAEEVVPNECSVHEATPPSGNGPQKKSVDRSIQTVVSCLF NPENRLRNSVVTQDDYFKDKRVHHFRRSRAMLKSV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001190811
ORF Size:	945 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001190811.1</u> , <u>NP_001177740.1</u>
RefSeq ORF:	948 bp
Locus ID:	1618
UniProt ID:	<u>Q92904</u>
Cytogenetics:	3p24.3
MW:	35.7 kDa