

Product datasheet for **RC231169**

Clorf201 (STPG1) (NM_001199012) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Clorf201
Synonyms:	Clorf201; MAPO2
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC231169 representing NM_001199012
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACAACCTCTGCACAGAAAAATGAACGCACTGGCAAACATCCCAGACGTGCCAGTGAAGTACAGAAAG
GTTTTACTGCTGCATATCCAACACAATCCTCCATTCTTTTAAATCTCAAGCTTCAGTAATCCCAGAATC
AGAAAAAAGGATTCAATAGTCAAGCCAAGAGATTTCTCATAAGAAGAATGATATCCCAGGACCTGGG
TTCTACAATGTTATTCACCAAGTCACCGGTGTCCAACAGTGTCTATTGTCCAAGAAAGGAAGTTCATGT
TTCCCTCAATGTGCGCCCGATTGGACACCATCATTTCTAAATACCCTGCAGCGAATGCATACACTATCCC
ATCGGATTTTATTTCCAAGAGAGACTTTAGTAATTCGTGTTCCAGCATGTTCCAGTTGCCAAGCTTTATG
AAAGCTCTCAAGTTTGAACTCCTGCACCAAACTATTACAATGCCTCTGTCTTGTCTGCAAGCAGAGAA
ACAACGTCTGTACTCGAGCCGGTTTATGTCAAAAACCCAAAGAGGATCTTTCGCTTTTGTGATAAAGG
ACCTCCCCCAGGGCATTATGATATCAACGAATCCCTTGTGAAGCAGTCGCCAAATACATTAATGTCTTGT
TTTAAATCAAAAACCAACCGTGATTTAACTGACGTCAACAGGCCCGGGACCTGGTTATTACAACCCCA
GTGATTGCACAAAAGTTCCAAAAAGACTCTTTCCCGAAAAACCCCATCTGAACTTCTCTGCTCAGCC
TTCGCTCTGCCTCCGAAGCCACCTTTCCAGGTCCTGGTCAGTATGAGATCGTGGACTACTTAGGCCCC
CGCAAGCATTTCATCTAGTGCATCATTCGTGTCCAATACCAGCCGGTGGACAGCGCGCCGCTCAGC
CAGGCTGCCTGGCCAGCTACGTACAAGCCAGAGCTTCCAGGAAAGCAGTCCTTCTCTACAACGAGGA
CAAGAAATGGATCCCGTTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATACGACGATAAGGTTTAA



Protein Sequence: >RC231169 representing NM_001199012
Red=Cloning site Green=Tags(s)


MDNSAQKNERTGKHPRRASEVQKGFTAAYPTQSSIPFKSQASVPISEKKGFNSQAKRFPHKKNDIPGPG
FYNVIHQSPVSNVSLSKKGTCTMFPSMCLRLDTIISKYPAANAYTIPSDFIISKRDFSNSCSCSMFQLPSFM
KALKFETAPAPNYNASVCCYQQRNNVCTRAGFMSKTRGQSFAFADKGGPPPHGQVINESL VKQSPNTLMSC
FASKTNRGLKLTSTGPGVGYNNPDSCTVKPKKTLFPKNPILNFSAPQSPHLPKPPFPGQYIEVDYLG
RKHFISASFSVNTSRWTAAPPOLGPGATPYKPELPGKGSFLYEDNKWKIPLV

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8052_h08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001199012

ORF Size: 1002 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001199012.1](#), [NP_001185941.1](#)

RefSeq ORF: 1005 bp

Locus ID: 90529


UniProt ID: [Q5TH74](#)

Cytogenetics: 1p36.11


MW: 37.2 kDa

Gene Summary: May positively contribute to the induction of apoptosis triggered by O(6)-methylguanine. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC231169



Western blot validation of overexpression lysate (Cat# [LY434168]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC231169 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).