

Product datasheet for **RC231180**

TEDC1 (NM_001198983) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	TEDC1
Synonyms:	C14orf80
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC231180 representing NM_001198983
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCGGCCAGCGCGGTGATTGGGTACAGACCTCCGCGCTCTGGCAGCTCCTCTTCGGTGTGCTCTCGC
CACTCCCTGCGGGCAACGCCTTGGCATCGCTCGCCCTGGAGGTCCAAGCCCGCTTGGTGAAGTCAGCACT
ATGCTCCAGGGCTACCCGAGGCTGGCACTGGCACAACCTACCTGAGGATGGCTCGCAGGGCAGTCGGGAG
CTGCTGCTGGCTCTGTCTGGCTCTTGGCCCGAGGACCTGTGCCCGAGCAGATGCTGGCCAGGCCCGAG
TGCCCTCTGGGTGACGAGATGACTGTGTGCCAGTGTGAGGCCCTGGCCAGCCCTGGCCACCTGCACCCCA
CATGGAAGCAGAGGGTCTGTGGATGTCCGCCATGTGCACTGGCTGATGGGAAAGCTGCGGTTCCGGTGG
CGCCAGCTGGTGTCCAGTCAGCAGGAGCAGTGCGCCCTCCTGAGCAAGCTGCTGCGGACTCTGGAGCGTG
AGAACCAGCGCCTGGAGGCTGTCTGGCGTGGCGGCGCTCTGAGCTGGTCTTCTGGCGGTGGATGGACAC
GGTCTCTGGGCACCTGTGCCCCGGAGGTGCCTGTGTCAGCCTCACAGCCACCTTCTGCCCTGGGTCCCC
GAGCGCGGGGGTGGCGAGTTGGACCTGGTAGTGCGGGAGCTGCAGGCACTGGAGGAGGAGCTGCGGGAGG
CTGCGGAGCGCAGGCGGGCGGCTGGGAGGCCAAGGCTGGAGGCTGTGGACGGGGGCCAGAGTGGAGTGC
CGCGCGGGCGGCTCTCGGGAGGCTGTGAAAAGGAGCTGGGAGCTCTACAGCAGTGTCTGGAGCGAGAC
GGTGGCCCGGCCAGCCCATGGGCCACACCGGCTGGTGAGACGAGAGGATGGGCGAGCAGGGGACCGGG
ACCTGCGGGCAGCTGTGGTGATCAGGACGCTGAGGAGCCAGGAGGCTGCCTGGAGGCGGTGCTACGTCTG
ACTACAGGACAGTGTGCGCAGGAAGTGGCCAGGCTGGTGGGAGCCGCCCTGGTCTCATCTGGATCCCC
CCACCTGGACGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC231180 representing NM_001198983

Red=Cloning site Green=Tags(s)

MSGQRGDWVQTSALWQLLFRVLSPLPAGNALASLAEVQARLVKSALCSQGYPRLLAQLPEDGSQGSRE
LLLALSWLLARGPVPEQMLAQARVPLGDGMTVCQCEALASPGPPAPHMEAEGPVDVRHVQWLMGKLRFRW
RQLVSSQEQCALLSKLLRTLRENQRLEAVLAWRRSELVFWRWMDTVLGTCAPEVPAASQPTFLPWVP
ERGGGELDLVVRELQALEEELREAAERRRAWEAKAGGCGRGPWSAARRASREAVEKELGALQQCWERD
GGPAQPHGPHRLVRREDGAAGDRDLRAAVVIRTLRSQEACLEAVLRRLLQGQCRQELARLVGARPGLIWIW
PPGR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8054_d12.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001198983

ORF Size:

1062 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:


The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001198983.1, NP_001185912.1</u>
RefSeq ORF:	1065 bp
Locus ID:	283643
UniProt ID:	<u>Q86SX3</u>
Cytogenetics:	14q32.33
MW:	39.7 kDa
Gene Summary:	May play a role in counteracting perturbation of actin filaments, such as after treatment with the actin depolymerizing microbial metabolite Chivosazole F.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC231180



Western blot validation of overexpression lysate (Cat# [LY434179]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC231180 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).