

## Product datasheet for **RC231181L3V**

### Aprataxin (APTX) (NM\_001195248) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Aprataxin (APTX) (NM_001195248) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Aprataxin
Synonyms:	AOA; AOA1; AXA1; EAOH; EOAHA; FHA-HIT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001195248
ORF Size:	1068 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231181).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001195248.1</a> , <a href="#">NP_001182177.1</a>
RefSeq ORF:	1029 bp
Locus ID:	54840
UniProt ID:	<a href="#">Q7Z2E3</a>
Cytogenetics:	9p21.1
Protein Families:	Druggable Genome
MW:	41.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the histidine triad (HIT) superfamily. The encoded protein may play a role in single-stranded DNA repair through its nucleotide-binding activity and its diadenosine polyphosphate hydrolase activity. Mutations in this gene have been associated with ataxia-ocular apraxia. Alternatively spliced transcript variants have been identified for this gene.[provided by RefSeq, Aug 2010]