

## Product datasheet for **RC231184L2V**

### GPR120 (FFAR4) (NM\_001195755) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR120 (FFAR4) (NM_001195755) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR120
Synonyms:	BMIQ10; GPR120; GPR129; GT01; O3FAR1; PGR4
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001195755
ORF Size:	1083 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231184).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001195755.1</a>
RefSeq ORF:	1086 bp
Locus ID:	338557
UniProt ID:	<a href="#">Q5NUL3</a>
Cytogenetics:	10q23.33
Protein Families:	Druggable Genome, Transmembrane
MW:	40.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a G protein-coupled receptor (GPR) which belongs to the rhodopsin family of GPRs. The encoded protein functions as a receptor for free fatty acids, including omega-3, and participates in suppressing anti-inflammatory responses and insulin sensitizing. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]