

Product datasheet for **RC231187L3V**

C10orf81 (PLEKHS1) (NM_001193434) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C10orf81
Synonyms:	C10orf81; HEL185
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001193434
ORF Size:	1098 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC231187).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001193434.1
RefSeq ORF:	1101 bp
Locus ID:	79949
UniProt ID:	Q5SXH7
Cytogenetics:	10q25.3
MW:	41.4 kDa

