

Product datasheet for **RC231205L3V**

C15orf38 (ARPIN-AP3S2) (NM_001199058) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ARPIN-AP3S2
Synonyms:	ARPIN; C15orf38; C15orf38-AP3S2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001199058
ORF Size:	1182 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC231205).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001199058.1
RefSeq ORF:	1185 bp
Locus ID:	100526783
Cytogenetics:	15q26.1
MW:	44.3 kDa



Gene Summary:

This locus represents naturally occurring read-through transcription between the neighboring C15orf38 (chromosome 15 open reading frame 38) and AP3S2 (adaptor-related protein complex 3, sigma 2 subunit) genes. The read-through transcript encodes a fusion protein that shares sequence identity with each individual gene product. [provided by RefSeq, Nov 2010]