

Product datasheet for RC231246

LOC646508 (NM_001195187) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	LOC646508
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC231246 representing NM_001195187
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCCGTCATTCCGTCCATCATCAGGCTCAGAGGCCCCAGAGCCAAAACCCACAGGAGCAGCAGA
GGAGACCCATGGGCGAGCAGCTCTCCAGCAGAGCAGGAGGAATCCAGGGTGAATGCAAGAACTGCGG
GGCCTTTGGGCACTCAGCCAGGAGCAAGACCTGCCCATTAAGAGGTGGAGTGGGGCCCTTCTCTGCAG
GCCCTGGGCTCACACAAGGAGAAGGAGAACCTGAAACCAGCAAGGCCAGCTACCCCTTTACGACTCCAG
GGCCCTTTACGACGAATGACAGAGAAAAGGAGCGAAGTCCAAGTCCCCAGCAGCAGCAGAGCGAAGCTCC
GACGCAGACATTTCCAGAACTCCCAAGAGAAAATGCAGGAAGCCTGGAAGGAGCCAGCAGAAGATTGT
TTGTTCTGAGGCATCCTACCATGCCACTGCCTGTCCACACCACCAAGAAGAGATCTGTCTGGGCCCTG
TGTCACAGGTCCACCGCTGTCAACAAACCCGAGATGAGATTACTCTGCCCTTCGGGTCAACACGATT
ACCTCAACTGAGCACCTGTGGACCCACCAAGGACATGGCAGGGACGTTACTGCCTCCCTGCTCCTGTT
CTGAAGAGTCCCACAGACCCCACTCTCAGTGCCAGGCTGCCAGCCAACAGGCTGACATGTCCTCCC
ATGGTGCTCTCCAGCCTGCCATGCAGGCGCTTGCCTGGGTCTGGCCTTAAATCCCAGGCAGAAATCAA
ACATCCCGACGCAGATGCAAAGCCCAGACCACAGCAAGTCAGAAAACAGTGTGGCCAGGACTCCAGAACC
CAGGCACAGACAAGGAGCCTGCCCCCGTCCCCACCCAGACTTTCAGAACCCCGCAAGAAAGCAAGAT
TCAGCTCCTTCCAGACCCCTGCACTGAGAACTCAGCTCCCGGATGTGGGCGCTGTGCAGACACTCCAGCC
TCCCCGCACTGCAACTGGACTTGGATCCAAAGAGGCACCCAAAGGCGACCCGAGAGACAGCAGCCACCAAG
ACAGCAACCTGCAAGCCAGAGTCAACCTCCAGCCGACCCAGCTCACCTTCTGGGCCAGCCAGG
GCTGCCCGCTCCTCCAGCTGGACACCCATTATGTTCCAGGGAGGCCGGCAGTGTACCTTCATGAG
AGGGGACAAGGGACAGAAGAGCCCCAGGTTCAGAATGCCTCCACATCCGTCCTCCTGAAAACCTGCT
TCTGCTCAGAGCCCTCGCTTCTCAAGGCAGCCTGAGGGGAGGGTCCCCAGGTCTCAACGAGTGTCTCT
ATGAGGACCTTCTGGTCACTTCTCTCTGAGGACAGTGACAGTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC231246 representing NM_001195187

Red=Cloning site Green=Tags(s)

MARHSVHHQAQRPPRAKNPQEQQRRPMGTTLPAEQEESRVKCKNCGAFGHSARSKTCPIKRWSGALPLQ
ALGSHKEKENLKPAKAQLPFTTPGPFTTNDREKERSPPQQQSEAPTQTFPRTPQEKMQEAWKEPAEDC
LFLRHPTMPLPVHTTKKRSVLGPVSTGPPPVNKPEMRLLCPSGHNDSPQLSTCGPTKGHGRDVTASLLPV
LKSSHQTPTLSARLPANRPDMSSHGALQPAMQALALGPGLKSQAEIKHPDADAKPRPQQVRKQCGQDSRT
QAPDKEPAPVPTQTFQNPAAKKARFSSFTPALRTQLPDVGAVQTLQPPRTATGLGSKEAPKATAETAATK
TATLQPRVNLQPAPSSPFLGPAQGGCPVLQPGPIHVPGRPGSVTFMRGDKGQKSPFRMPPTSRRPENSA
SAQSPRFSRQPEGQGPQVSTSVLYEDLLVTSSSDSDSD

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/ja1306_g04.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001195187

ORF Size:

1377 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001195187.1](#), [NP_001182116.1](#)


RefSeq ORF: 1379 bp

Locus ID: 646508


Cytogenetics: 19q13.42

MW: 50.1 kDa

Product images:



Circular map for RC231246



Western blot validation of overexpression lysate (Cat# [LY434245]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC231246 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).