

Product datasheet for **RC231286L3V**

eIF4GII (EIF4G3) (NM_001198803) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | eIF4GII (EIF4G3) (NM_001198803) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | EIF4G3 |
| Synonyms: | eIF-4G 3; eIF4G 3; eIF4GII |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001198803 |
| ORF Size: | 1545 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC231286). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001198803.1 |
| RefSeq ORF: | 1548 bp |
| Locus ID: | 8672 |
| UniProt ID: | O43432 |
| Cytogenetics: | 1p36.12 |
| Protein Families: | Transcription Factors |
| Protein Pathways: | Viral myocarditis |
| MW: | 55.7 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is thought to be part of the eIF4F protein complex, which is involved in mRNA cap recognition and transport of mRNAs to the ribosome. Interestingly, a microRNA (miR-520c-3p) has been found that negatively regulates synthesis of the encoded protein, and this leads to a global decrease in protein translation and cell proliferation. Therefore, this protein is a key component of the anti-tumor activity of miR-520c-3p. [provided by RefSeq, May 2016]