

## Product datasheet for **RC231322L4V**

### Myoneurin (MYNN) (NM\_001185119) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Myoneurin (MYNN) (NM_001185119) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Myoneurin
Synonyms:	OSZF; SBBIZ1; ZBTB31; ZNF902
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001185119
ORF Size:	1743 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231322).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001185119.1</a> , <a href="#">NP_001172048.1</a>
RefSeq ORF:	1746 bp
Locus ID:	55892
UniProt ID:	<a href="#">Q9NPC7</a>
Cytogenetics:	3q26.2
Protein Families:	Transcription Factors
MW:	65.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the BTB/POZ and zinc finger domain-containing protein family that are involved in the control of gene expression. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 14. [provided by RefSeq, Jun 2010]